**ARGON MEDICAL DEVICES, INC.**  
 1445 Flat Creek Road  
 Athens, Texas 75751 USA  
 Tel: +1 800-927-4669  
 Tel: +1 903-675-9321  
 www.argonmedical.com

**EC REP EMERGO EUROPE**  
 Westervoortsedijk 60  
 6827 AT Arnhem  
 The Netherlands  
 +30 70 345 8570

**UK REP Emergo Consulting (UK) Limited**  
 c/o Cr 360 - UL International  
 Compass House, Vision Park Histon  
 Cambridge CB24 9BZ  
 United Kingdom

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Brazilian Portuguese	O glossário dos símbolos está localizado eletronicamente em <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Bulgarian	Речникът на символите е наличен в електронен вид на адрес <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Czech	Slovník symbolů v elektronické formě je umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Danish	Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
German	Das Glossar der Symbole kann unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Spanish	El glosario de símbolos puede consultarse en línea en la dirección <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Finnish	Elektroninen symbolisanasto on osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
French	Le glossaire des symboles est disponible au format électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Lithuanian	Simbolių žodyno elektroninė forma pateikta <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Norwegian	Symbolforklaringen befinner seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Romanian	Glosarul simbolurilor este stocat electronic la <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Swedish	Symbolforklaringen hittas elektroniskt på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>

Turkish	Semboller sözlüğü elektronik ortamda <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> adresinde yer almaktadır
Arabic	يتم توفير مصدر الرموز إلكترونيًا على <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>

IFU1100-0221  
 Date of Issue: 08/18/2023 Rev. B

en - English

**WORKER™ Guidewire, Amplatz and Bentson Type**

**Intended Use/Purpose**

The device is used to facilitate the placement of catheters.

**Device Description**

WORKER™ guidewires are made of stainless steel core with a polytetrafluoroethylene (PTFE) coated stainless steel spring wire that helps reduce friction during procedure. It offer a safe and flexible guidewire.

**Indications for use**

WORKER™ Guidewires are for use in facilitating the placement of catheters within the coronary and peripheral vasculature. They are also for use in the placement of non-vascular catheters.

**Duration**

Transient, less than 60 minutes.

**Warnings**

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package appears open or if the expiry date has been exceeded.
- Do not advance the wire against resistance until the cause of the resistance has been determined by fluoroscopy. Excess force against resistance may result in damage to guidewire, or catheter, or vessel perforation.
- Do not excessively torque the guidewire.
- Do not entrap or over-rotate the distal tip of the guidewire, it can lead to wire rupture.
- Do not excessively bend the guidewire, it may result in wire fracture.
- Guidewire placement should be monitored during placement and manipulation using fluoroscopic or suitable imaging method.
- Do not withdraw a guidewire through a needle. Straighten the guidewire in order to withdraw the needle.
- The product must only be used by qualified personnel who are familiar with the technique.

**Precaution**

- Twisting or entanglement of the guidewires may occur when more than one wire is being used simultaneously. This can be prevented by carefully isolating and marking the proximal end of the wires.

**Potential Complications**

Potential complications include, but are not limited to the following:

- Perforation of vessel or arterial wall
- Thrombus formation
- Puncture site hematoma
- Vasospasm
- Infection
- Vessel dissection

**Preparation**

- Make sure that the guidewire is undamaged.
- Loosen the straightener together with the guidewire.
- Inject saline with syringe in a vertical position.
- Make sure the product is compatible with other equipment.

**Directions for Use**

- Insert the flexible end of the guidewire into the device. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device.
  - If "J" wire is used, prepare "J" tip by sliding J-straightener over tip to straighten. Insert J-tip inside J-Straightener into device hub gently.
- Advance the guidewire through the device.
- Position the guidewire.
- Hold the guidewire in place and withdraw the device.
- Advance the catheter over the guidewire.

**Storage**

Store at controlled room temperature. Avoid extended exposure to light.

**Disposal**

After use, handle and dispose in accordance with end-user facility's policies and procedures concerning biohazardous materials and waste.

**NOTE :** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at [quality.regulatory@argonmedical.com](mailto:quality.regulatory@argonmedical.com) as well as to the competent health authority where the user/patient resides.

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