


Quick Reference Guide in GI

Primary Application Guidance


Oozing bleed or resected area ⁽¹⁾



**STEP
01**

Clear


Remove as much blood and fluid as possible



**STEP
02**

Direct


Apply PuraStat® at the base of the lesion, as close as possible to the source of bleeding



**STEP
03**

Apply

Apply PuraStat in an even layer until the whole lesion is covered



Tip: Note the direction of the flow (gravity), start at the uppermost point and work downwards.

**STEP
04**

Wait

Do not disturb PuraStat until sufficient time is allowed for haemostasis

Re-Application Guidance

If necessary, repeat the application

**STEP
05**

Direct

Place applicator through the previous PuraStat layer directly onto source of bleeding

**STEP
06**

Reapply

Apply PuraStat at the base of the lesion, as close as possible to the source of bleeding

**STEP
07**

Wait

Do not disturb PuraStat until sufficient time is allowed for haemostasis

**After
haemostasis
is confirmed,
excess PuraStat
may be left in place**

Indications

PuraStat® is indicated for haemostasis in the following situations encountered during surgery, when haemostasis by ligation or standard means is insufficient or impractical:

- Bleeding from small blood vessels and oozing from capillaries of the GI tract following surgical procedures.
- PuraStat is also indicated for the reduction of delayed bleeding following gastrointestinal endoscopic submucosal dissection (ESD) procedures in the colon.



3-D Matrix Europe SAS

11 chemin des Petites Brosses
69300 Caluire-et-Cuire-FRANCE
Tel: +33 (0)4 27 19 03 40
infoeu@puramatrix.com
www.3dmatrix.com


Reference:
(1) PuraStat IFU-002 Rev 2.2

PuraStat is a class III medical device, CE marked according to European council directive 93/42/EEC on medical devices and its relatives



Catheter

- Attach the syringe to catheter using the **luer lock connectors**.
- Ensure the connection is **finger tight**.
- Push the plunger on the syringe to **slowly prime** the catheter.
- Stop as soon as PuraStat reaches the end of the catheter.
- It is now **ready for use**.
- Estimated priming volumes
 - 2200mm Catheter Type-E (REF 634142): 2.2ml
- To utilize all of PuraStat in the **catheter's dead space**, remove the syringe, draw up room air, reattach and slowly push through the remaining PuraStat.



TOP Corporation
19-10 Senjunakai-cho
Adachi-ku, Tokyo
120-0035 Japan

PuraStat Nozzle System Type-E, is a class IIa medical device, CE marked according to European council directive 93/42/EEC on medical devices and its relatives



PuraStat Nozzle System type E is distributed by
3-D Matrix Europe SAS

11 chemin des Petites Brosses
69300 Caluire-et-Cuire-FRANCE
Tel: +33 (0)4 27 19 03 40
infoeu@puramatrix.com
www.3dmatrix.com



DUOMED

Duomed Swiss AG
Grenzstrasse 5a, 6214 Schenkon
Tel. +41 (0)41 510 07 00
www.duomed.com

EC REP MPS Medical Product Service GmbH

Borngasse 20
35619 Braunfels
Germany