

# Disposable 3-Ply 3-folder ear-loop face masks with shield



REPEL mask is a uniquely designed pleated 3-ply, 3 folder tie-on face mask, with anti fog shield or fog free foam strip to provide extra protection. It is 100% latex-free and fiberglass free. All materials are non-irritant and odorless. Available in both tie-on and ear-loop type.

## With Earloops/Ties

- Complies with EN14683 Type IIR
- BPF  $\geq 98\%$  at 3 microns
- PFE  $\geq 98\%$  at 0.1 micron
- Resistance to penetration by synthetic blood - 160mm/Hg

CE  LATEX FREE NICKEL FREE CADMIUM FREE

- 3-ply, 3-folder comfortable design.
- Fluid resistant perforated film (PPF) outer layer.
- Meltblown filtration middle layer.
- Fluid resistant perforated film (PPF) as third layer can provide additional protection and breathability.
- Silk-like, ultra-soft cellulosic white inner layer.
- 100% latex, fiberglass free hypoallergenic material.
- Full length aluminum nose piece.

Item No.	Description	Packing	Size of Case
BK1400D	3-Ply 3-folder ear-loop face masks with shield	25pcs/box, 10boxes/case	40.4*33.4*24 CM

**A:** FLUID RESISTANT PERFORATED FILM (PPF)

**B:** : SILK LIKE, ULTRA SOFT INNER LAYER

**C:** FLUID RESISTANT PERFORATED FILM (PPF)

