



DUOMED

# Endoscopy Catalogue

**Bringing high quality**

products, services

& innovations

to the healthcare

community



# Table of contents

## Duomed UK

Introduction	4
Training Centre & Skills Lab	5

## Endo Essentials

### Protection & Cleaning

Water Jet Connector	7
Biopsy Valves	8
Air/Water Valves	9
Suction Valves	10
Cleaning Valves	11
Biopsy Valve Sets	12
Bite Blocks	13

### Tissue Acquisition

Standard Biopsy Forceps	14
Alligator Biopsy Forceps	15
Jumbo Biopsy Forceps	16
Hot Biopsy Forceps	17
Tracheo-Bronchial Biopsy Forceps	18
Paediatric Biopsy Forceps	19
Cytology Brushes	20
Pulmonary Brushes	20

### Patient Clothing

Dignity Shorts	21
----------------	----

### Cleaning & Disinfection

EcoFlushCaps with CO2	22
EcoFlushCaps without CO2	23
HybriX EcoFlushCaps with CO2	24
Irrigation Tubing	25
CO2 Tubing	26
Jet Channel Adapters	27
Lubricating Jelly	28

## Endo Therapy

### Diagnostic Colouring

Lugol	30
Acetic Acid	30
Indigo Carmine	30
Spray Catheters	31

### Dilation

Dilation Balloons	32
Dilation Balloon 3 Diameters	33
Disposable Bougies	34
3 Stage Wire Guided Balloons	35

Inflation Devices	36
-------------------	----

### EMR

Endoscopic Hood	37
DiLumen EZ Glide	38

### ESD

Splash M-Knife	39
ESD Knives	40
Bipolar HemoStat WideCup	41
DiLumen Dynamic Retraction	42

### ERCP

Sphincterotomes	43
Stone Extraction Baskets	44
Stone Retrieval Balloons	45
Surgical Drape	46

### Haemostasis

NovaClip R2	47
NovaClip R4	48
HemoClips	49

### Injection Needles

Novaject EMR Needles	50
Injection Needles	51

### Snares

Cold Snares	52
Rotatable Snares	53
Monofilament Snares	54
Multi-Stage Snares	55
Oval Hot Snares	56
Hexagonal Hot Snares	57
Crescent Hot Snares	58
Round Snares	59
Braided Snares	60
Multi-Opening Snares	61
K-Snare	62

### Retrieval Management

Retrieval Nets	63
Polyp Traps	64
Grasping Forceps	65

## Endo Decontamination

### Cleaning & Disinfection

Single-Ended Cleaning Brushes	67
Double-Ended Cleaning Brushes	68
Multi-Point Cleaning Brushes	69
Triple-Point Cleaning Brushes	70

Sweeper Brushes	71
-----------------	----

Combination Cleaning Brushes	72
------------------------------	----

X-MAS Cleaning Brushes	73
------------------------	----

Valve Brushes	74
---------------	----

U-Brush	75
---------	----

Cleaning Brush Sets	76
---------------------	----

### Transport

Transport Bags	77
----------------	----

### Protecion & Cleaning

Distal Protectors	78
-------------------	----

3-in-1 Cleaning Sponges	79
-------------------------	----

## GI Stenting

Oesophageal Stents	81
--------------------	----

Colonic Stents	82
----------------	----

Standard Guidewires	83
---------------------	----

Special Guidewires	84
--------------------	----

## Capsule Endoscopy

Small Bowel Capsule Endoscopy	86
-------------------------------	----

Stomach Capsule Endoscopy	87
---------------------------	----

## Medi Tech

POSE 2.0	89
----------	----

EndoVE	90
--------	----

CUTIS	91
-------	----

CADDIE	92
--------	----

## Clinical Training

The TD Box	94
------------	----

## Watchman Managed Inventory

	96
--	----

# Introduction



Duomed UK is an innovative organisation, breaking new ground in the endoscopy healthcare environment. We are driven to provide outstanding customer service and product solutions.

All our teams are united in the understanding that what we do every day, makes a difference to our customers and their patients.

We apply purpose driven products and solutions from around the globe, to help facilitate improved delivery of healthcare in the NHS and private sectors.

Duomed UK is a part of The Duomed Group, a dynamic organization with a well-established reputation. The Duomed Group are active in consultancy, sales, integration, training and technical support of medical technology and devices for hospitals and medical practices.

Our group currently consists of Acertys, Braun Scandinavia, Duomed, Dyma Medical, FMH Medical, Duomed UK, Life Partners Europe and Treier Endoscopie.

Discover all our products, services, and courses at [www.duomed.com/uk](http://www.duomed.com/uk).



# Clinical Training Suite

The Duomed UK Head Office has its own Clinical Training Suite, equipped with portable facilities and a wide range of accessories for in-house virtual, and hands on training.

Our Clinical Training Suite has multiple workstations, allowing groups to train using our T-D Training Boxes and utilising our extensive product portfolio.

While animal models are used for advanced training, Duomed UK also offer the option of a purpose-designed manufactured silicone model, perfect for the early stages of training or for staff who may have religious, ethical or personal objections to using animal models.

For any questions or enquiries please contact us at [info.uk@duomed.com](mailto:info.uk@duomed.com).

We look forward to welcoming you to our Clinical Training Suite.



# ENDO ESSENTIALS



# Water Jet Connectors

## Protection & Cleaning



- Minimizes the risk of cross contamination
- Eliminates the need for manual reprocessing

Supplier	Product Code	Description	Compatible Endoscope	QOM	FAL Code
Vytil	VTL/CN03	Water jet connector	Fujifilm	100	
Vytil	VTL/CN02	Water jet connector	Pentax	100	

# Biopsy Valves



CJ-VAB-OF



NLS/BVNS-2-OFP



BV-PX-01



BV-OX-01



BV-FX-02-EO

- Minimizes the risk of cross contamination
- Eliminates the need for manual reprocessing

Supplier	Product Code	Description	Compatible Endoscope	QOM	FAL Code
Vytil	CJ-VAB-OF	Biopsy Valve Diam. 6.3mm Height 11.3mm	Olympus / Fujifilm	50	FAL24657
Vytil	CJ-VAB-P	Biopsy Valve Diam. 7.3mm Height 17.3mm	Pentax	50	FAL24648
Vytil	NLS/BVNS-2-OFP	3-in-1 Biopsy Valve	Olympus / Fujifilm / Pentax	50	
Vytil	VTL/BP01-20	Biopsy Valve Non Sterile Single Use Pouches of 20	Olympus / Fujifilm	600	
Vytil	VTL/BP01S	Biopsy Valve Single Use	Olympus / Fujifilm	100	
Vytil	VTL/BP02S	Biopsy Valve Single Use	Pentax	100	
Vytil	VTL/BP0701S	Biopsy Valve Irrigator Single Use	Pentax	100	
Vytil	VTL/BP07S	Biopsy Valve Irrigator Single Use	Olympus / Fujifilm	100	
MDD	BV-FX-01	Biopsy Valve Non Sterile Single Use Green	Fujifilm	200	FAL23563
MDD	BV-FX-02-EO	Biopsy Valve Sterile Single Use Green	Fujifilm	200	FAL23618
MDD	BV-OX-01	Biopsy Valve Non Sterile Single Use Blue	Olympus	200	FAL23723
MDD	BV-OX-02-EO	Biopsy Valve Sterile Single Use Blue	Olympus	200	FAL23520
MDD	BV-PX-01	Biopsy Valve Non Sterile Single Use Red	Pentax	200	FAL23639
MDD	BV-PX-02-EO	Biopsy Valve Sterile Single Use Red	Pentax	200	FAL23737

Water Jet Connectors

Biopsy Valves

# Air/Water Valves



AW-FJ-01



AW-OLY-01

- Minimizes the risk of cross contamination
- Supports constant performance

# Suction Valves



SU-FJ-01



SU-OLY-01

- Minimizes the risk of cross contamination
- Standardises infection prevention practices in the GI unit
- Supports constant performance

Supplier	Product Code	Description	Compatible Endoscope	QOM	FAL Code
MDD	AW-FJ-01	Air/Water Valve	Fujifilm	100	
MDD	AW-OLY-01	Air/Water Valve	Olympus	100	

Supplier	Product Code	Description	Compatible Endoscope	QOM	FAL Code
MDD	SU-FJ-01	Suction Valve	Fujifilm 700	100	
MDD	SU-OLY-01	Suction Valve	Olympus	100	



AWV-KIT-07-FJ



AWV-KIT-07-OLY

- Minimises risk of cross-contamination
- Supports constant performance

Supplier	Product Code	Description	Compatible Endoscope	QOM	FAL Code
MDD	AWV-KIT-07-OLY	Cleaning Valve	Olympus	100	FAL23568
MDD	AWV-KIT-07-FJ	Cleaning Valve	Fujifilm 700	10	FAL23815



OLYMPUS SET



FUJIFILM SET

- Minimises risk of cross-contamination
- Single use & individually sterile packed

Supplier	Product Code	Description	Compatible Endoscope	QOM	FAL Code
MDD	AWV-KIT-01-OLY	Biopsy Air/Water Suction, Water Jet Adapter	Olympus	25	FAL23578
MDD	AWV-KIT-02-OLY	Biopsy Valve, Air/Water Valve, Suction Valve	Olympus	25	FAL23708
MDD	AWV-KIT-03-OLY	Air/Water Valve, Suction Valve	Olympus	25	FAL23722
MDD	AWV-KIT-04-OLY	Biopsy Air/Water Suction, Air/Water Cleaning	Olympus	25	FAL23749
MDD	AWV-KIT-05-OLY	Biopsy A/W Suction A/W Cleaning + Jet Adapter	Olympus	25	FAL23818
MDD	AWV-KIT-06-OLY	Air/Water Suction + Air/Water Cleaning	Olympus	25	FAL23908
VYTIL	NLS/KBVNS-O	Biopsy Valve, Suction Valve, Air/Water Valve	Olympus	30	FAL23587
MDD	AWV-KIT-01-FJ	Biopsy Valve, Air/Water Valve, Suction Valve	Fujifilm 700	25	FAL23664
MDD	AWV-KIT-02-FJ	Air/Water Valve, Suction Valve	Fujifilm 700	25	FAL23558
MDD	AWV-KIT-03-FJ	Biopsy Valve, Air/Water Valve, Suction Valve	Fujifilm 500/600	25	FAL24602
MDD	AWV-KIT-04-FJ	Air/Water Valve, Suction Valve	Fujifilm 500/600	25	FAL24602
MDD	AWV-KIT-05-FJ	Biopsy Air/Water Suction, Jet Adapter	Fujifilm 700	25	FAL24564
MDD	AWV-KIT-06-FJ	Biopsy Air/Water Suction + Cleaning Valve	Fujifilm 700	25	FAL24564
MDD	AWV-KIT-01-PX	Biopsy Air/Water Suction, Jet Adapter	Pentax	25	FAL24585
MDD	AWV-KIT-02-PX	Biopsy Valve, Air/Water Valve, Suction Valve	Pentax	25	FAL24607
MDD	AWV-KIT-03-PX	Air/Water Valve, Suction Valve	Pentax	25	FAL24584
MDD	AWV-KIT-04-PV	Biopsy Air/Water Suction + Cleaning	Pentax	25	FAL24584

# Bite Blocks



- To protect the endoscope during a gastroenteroscopy
- Strap material: textile, rubber, spandex or no strap
- Various sizes available



Supplier	Product Code	Description	Strap	Opening diameter mm / fr	QOM	FAL Code
Vytil	CJ-JHY-KD-B1	Bite Block	Rubber	20 / 60	25	
Vytil	NLS/BBA	Bite Block	Spandex	20 / 60	100	FAL23867
Vytil	NLS/BBS-O2	Oxygen Bite Block	Spandex	17.8 x 26 / 53	100	FAL23617
Vytil	VTL/BT07	Pediatric Bite Block	Spandex	19 x 26 / 48	50	FAL24653
MDD	BB-01-85-50-R	Bite Block	Rubber	27 x 23	100	FAL23840
MDD	BB-03-30-33-R	Pediatric Bite Block	Rubber	19 x 16	100	FAL23763
MDD	BB-05-30-30-R	Bite Block	Rubber	25 x 22	100	FAL23640
MDD	BB-05-30-OXY-R	Oxygen Bite Block	Rubber	25 x 18	100	FAL23821
MDD	BB-01-85-50-T	Bite Block	Textile	27 x 23	100	FAL23634
MDD	BB-03-30-33-T	Pediatric Bite Block	Textile	19 x 16	100	FAL23496
MDD	BB-05-30-30-T	Bite Block	Textile	25 x 22	100	FAL23489
MDD	BB-05-30-OXY-T	Oxygen Bite Block	Textile	25 x 18	100	FBE24586
MDD	BB-01-85-50	Bite Block	No Strap	27 x 23	100	FAL23809
MDD	BB-03-30-33	Pediatric Bite Block	No Strap	19 x 16	100	FAL23682
MDD	BB-05-30-30	Bite Block	No Strap	25 x 22	100	FAL23711
MDD	BB-05-30-OXY	Oxygen Bite Block	No Strap	25 x 18	100	FAL23533

# Standard Biopsy Forceps

## Tissue Acquisition



Supplier	Product Code	Compatibility	Spike	Catheter		QOM	FAL Code
				Length (cm)	Diameter (mm)		
Vytil	CJ-PPT-18-110-A	Tracheo-bronchial	Yes	110	1.8	10	FAL23566
Vytil	CJ-PAT-23-160	Gastric	No	160	2.3	10	
Vytil	CJ-PAT-23-160-A	Gastric	Yes	160	2.3	10	
Vytil	CJ-PAT-23-230	Gastric and Colonic	No	230	2.3	10	
Vytil	CJ-PAT-23-230-A	Gastric and Colonic	Yes	230	2.3	10	
Vytil	NLS/BFOS-24-160	Gastric and Colonic	Yes	160	2.4	100	
Vytil	NLS/BFOS-24-230	Gastric and Colonic	Yes	230	2.4	100	
Vytil	NLS/BFO-24-160	Gastric and Colonic	No	160	2.4	100	
Vytil	NLS/BFO-24-230	Gastric and Colonic	No	230	2.4	100	
MDD	BF12-ECO-18-180-SU	Bronchoscope	Yes	180	1.8	10	
MDD	BF11-ECO-18-180-SU	Bronchoscope	No	180	1.8	10	
MDD	BF12-ECO-18-230-SU	Bronchoscope	Yes	230	1.8	10	
MDD	BF11-ECO-18-230-SU	Bronchoscope	No	230	1.8	10	
MDD	BF12-ECO-23-160-SU	Gastric and Colonic	Yes	160	2.3	10	FAL23912
MDD	BF11-ECO-23-160-SU	Gastric and Colonic	No	160	2.3	10	FAL23873
MDD	BF12-ECO-23-180-SU	Gastric and Colonic	Yes	180	2.3	10	FAL23628
MDD	BF11-ECO-23-180-SU	Gastric and Colonic	No	180	2.3	10	FAL23565
MDD	BF12-ECO-23-230-SU	Gastric and Colonic	Yes	230	2.3	10	FAL23674
MDD	BF11-ECO-23-230-SU	Gastric and Colonic	No	230	2.3	10	FAL23512
MDD	BF11-ECO-S-23-160-SU	Gastric and Colonic	No	160	2.3	10	
MDD	BF12-ECO-S-23-180-SU	Gastric and Colonic	Yes	180	2.3	10	
MDD	BF11-ECO-S-23-180-SU	Gastric and Colonic	No	180	2.3	10	
MDD	BF12-ECO-S-23-230-SU	Gastric and Colonic	Yes	230	2.3	10	
MDD	BF11-ECO-S-23-230-SU	Gastric and Colonic	No	230	2.3	10	
MDD	BF12-ECO-18-280-SU	Entroscope	Yes	280	1.8	10	
MDD	BF11-ECO-18-280-SU	Entroscope	No	280	1.8	10	
Pentax	FB-18-105-W-P-B1	≥Ø 2.0	Yes	105	2	10	
Pentax	FB-18-105-O-P-B1	≥Ø 2.0	No	105	2	10	
Pentax	FB-18-160-O-P-B1	≥Ø 2.0	No	160	2	10	
Pentax	FB-23-230-W-P-B1	≥Ø 2.8	Yes	230	2.8	10	
Pentax	FB-23-230-O-P-B1	≥Ø 2.8	No	230	2.8	10	




# Alligator Biopsy Forceps



Supplier	Product Code	Compatibility	Spike	Catheter		QOM	FAL Code
				Length (cm)	Diameter (mm)		
Vytil	CJ-PATC-23-160	Gastric	No	160	2.3	10	
Vytil	CJ-PATC-23-160-A	Gastric	Yes	160	2.3	10	
Vytil	CJ-PATC-23-230	Gastric and Colonic	No	230	2.3	10	
Vytil	CJ-PATC-23-230-A	Gastric and Colonic	Yes	230	2.3	10	FAL24643
Vytil	NLS/BFC-24-160	Gastric	No	160	2.4	100	
Vytil	NLS/BFCS-24-160	Gastric	Yes	160	2.4	100	
Vytil	NLS/BFC-24-230	Gastric and Colonic	No	230	2.4	100	
Vytil	NLS/BFCS-24-230	Gastric and Colonic	Yes	230	2.4	100	
MDD	BF11-ALI-18-180-SU	Bronchoscope	No	180	1.8	10	
MDD	BF12-ALI-18-180-SU	Bronchoscope	Yes	180	1.8	10	
MDD	BF11-ALI-18-230-SU	Bronchoscope	No	230	1.8	10	
MDD	BF12-ALI-18-230-SU	Bronchoscope	Yes	230	1.8	10	
MDD	BF11-ALI-23-180-SU	Gastric and Colonic	No	180	2.3	10	FAL23598
MDD	BF12-ALI-23-180-SU	Gastric and Colonic	Yes	180	2.3	10	FAL23896
MDD	BF11-ALI-23-230-SU	Gastric and Colonic	No	230	2.3	10	FAL23484
MDD	BF12-ALI-23-230-SU	Gastric and Colonic	Yes	230	2.3	10	FAL23527
MDD	BF11-ALI-18-280-SU	Entroscope	No	280	1.8	10	
MDD	BF12-ALI-18-280-SU	Entroscope	Yes	280	1.8	10	FAL23613
Pentax	FB-23-230-O-P-A1	≥ 2.8	No	230	2.8	10	
Pentax	FB-23-230-W-P-A1	≥ 2.8	Yes	230	2.8	10	
Pentax	FB-18-105-O-P-A1	≥ 2.0	No	105	1.8	10	

# Jumbo Biopsy Forceps



- Graduated catheter to easily identify the biopsy forceps progression within the operating channel
- All forceps are delivered with a sample retriever to reduce blood exposure risk

Supplier	Product Code	Description	Spike	Catheter		QOM	FAL Code
				Length (cm)	Diameter (mm)		
Vytil	AN/AMHBFJ-230	Jumbo	No	230	2.8	100	
MDD	BF12-ECO-30-230-SU	Colonic	Yes	230	3	10	FAL23471
MDD	BF11-ECO-30-230-SU	Colonic	No	230	3	10	FAL23557


# Hot Biopsy Forceps



- This forcep is used for small polyps removal
- Oval and fenestred jaws made of surgical stainless steel
- PTFE coated catheter
- Coagulation is achieved with opened or closed jaws

Supplier	Product Code	Forceps		Catheter		QOM	FAL Code
		Jaws	Spike	Length (cm)	Diameter (mm)		
Vytil	CJ-PC-23-230	= standard	No	230	2.3	10	FAL23688
Vytil	CJ-PCA-23-230	= standard	Yes	230	2.3	10	FAL23583
Vytil	NLS/HotBFO-24-230	= standard	No	230	2.4	10	

# Tracheo-Bronchial Biopsy Forceps



Alligator jaw with spike

Standard jaw without spike


- Graduated catheter to easily identify the biopsy forceps progression within the operating channel
- All forceps are delivered with a sample retriever to reduce blood exposure risk

Supplier	Product Code	Forceps			Catheter		QOM	FAL Code
		Jaws	Spike	PTFE coated	Length (cm)	Diameter (mm)		
Vytil	AN/AMHBFA	= standard	No	Yes	110	1.8	10	
Vytil	AN/AMHBFA-WC	= standard	No	No	110	1.8	10	
Vytil	AN/AMHBFC	🐊 alligator	No	Yes	110	1.8	10	
Vytil	AN/AMHBFC-WC	🐊 alligator	No	No	110	1.8	10	
Vytil	CJ-PPT-18-110-A	= standard	Yes	Yes	110	1.8	10	FAL24646



- Markers on the catheters to allow for fast pushing through the scope
- Tissue sampler to ease tissue removal from the forcep

Supplier	Product Code	Forceps		Catheter		QOM	FAL Code
		Jaws	Spike	Length (cm)	Diameter (mm)		
Vytil	NLS/BFO-18-160	= standard	No	160	1.8	100	
Vytil	CJ-PPT-18-160	= standard	No	160	1.8	20	FAL24634
Vytil	CJ-PPT-18-160-A	= standard	Yes	160	1.8	20	FAL24633



Cytology Brush

Pulmonary Brush

### Cytology Brush

- Soft brush for atraumatic collection
- Bronchial mucosa cells cytology brush

### Pulmonary Brush

- Pulmonary brush with a double catheter to prevent the contamination of the sample
- Catheter tip operculum
- Delivered with a sterile pair of scissors

Supplier	Product Code	Description	Catheter Diameter (mm)	Usable Length (cm)	Brush Diameter (mm)	Brush Length (cm)	QOM	FAL Code
Vytil	CJ-BC-18-110-A2	Cytology Brush	1.8	110	2	10	10	FAL23526
Vytil	CJ-BC-18-110-A3	Cytology Brush	1.8	110	3	10	10	FAL23532
Vytil	CJ-BM-17-105-S2	Pulmonary Brush incl. scissors	1.7	105	2	10	10	FAL23529
MDD	CYB-18-04-10-120	Cytology Brush Straight	1.8	120	4	10	10	FAL23822
MDD	CYB-23-04-10-120	Cytology Brush	2.3	230	4	10	10	FAL23686
MDD	CYB-23-04-10-160	Cytology Brush Straight	2.3	160	4	10	10	FAL23589
MDD	CYB-23-04-10-180	Cytology Brush Straight	2.3	180	4	10	10	FAL23695
MDD	CYB-23-04-10-230	Cytology Brush Straight	2.3	230	4	10	10	FAL23523
MDD	CYB-US-18-04-10-120	Cytology Brush U-Shape	1.8	120	4	10	10	FAL23894
MDD	CYB-US-23-04-10-120	Cytology Brush U-Shape	2.3	120	4	10	10	
MDD	CYB-US-23-04-10-160	Cytology Brush U-Shape	2.3	160	4	10	10	FAL23888
MDD	CYB-US-23-04-10-180	Cytology Brush U-Shape	2.3	180	4	10	10	FAL23726
MDD	CYB-US-23-04-10-230	Cytology Brush U-Shape	2.3	230	4	10	10	FAL23846

# Dignity Shorts

## Patient Clothing



### Greener Endoscopy Dignity Shorts

- Made from 100% recycled material & can be decomposed naturally
- 50% less plastic packaging
- Made from a non-woven fabric, utilizing a new generation of environmental protection
- Made with a non-toxic material
- Flap with Velcro fastenings allowing for cover to be kept in the open or closed position.

### Vytal Dignity Shorts

- Back opening for scope access without removing the shorts
- Velcro scratch discreet closure

Supplier	Product Code	Short Length (cm)	Waist Circumference (cm)	QOM	FAL Code
Greener Endoscopy	MDDS-L	65	120	100	FAL23525
Greener Endoscopy	MDDS-XL	65	130	100	FAL23750
Vytal	COLOPANTS-L-BD	55	120	100	
Vytal	COLOPANTS-XL-BD	60	140	100	

# EcoFlushCaps with CO2

## Cleaning & Disinfection

mdd



- Universal cap that fits all standard water bottles (Baxter/Fresenius)
- Available with standard Co2 connector hose or 1.1m connection hose options
- Available for Pentax, Fujifilm and Olympus scopes

Supplier	Product Code	Description	Compatibility	QOM	FAL Code
MDD	WB-FJ3-WC2	with CO2 connection	Fujifilm 500/600	10	FAL24636
MDD	WB-FJ4-WC2	with 1.1m CO2 connection	Fujifilm 500/600	15	FAL23727
MDD	WB-FJ5-WC2	with 1.1m CO2 connection	Fujifilm 700	15	FAL23913
MDD	WB-FJ6-WC2	with CO2 connection	Fujifilm 700	10	FAL24660
MDD	WB-FJ7-WC2	Twin EcoFlushCap	Fujifilm 700	15	FAL24569
MDD	WB-OL3-WC2	with CO2 connection	Olympus	10	FAL24629
MDD	WB-OL4-WC2	with 1.1m CO2 connection	Olympus	15	FAL23838
MDD	WB-PX5-WC2	with 1.1m CO2 connection	Pentax	15	FAL23692
MDD	WB-PX6-WC2	with CO2 connection	Pentax	10	

## EcoFlushCaps without CO2



WB-OL2-NC2

- Universal cap that fits all standard water bottles (Baxter/Fresenius)
- Available for Pentax, Fujifilm and Olympus scopes

Supplier	Product Code	Compatibility	QOM	FAL Code
MDD	WB-FJ2-NC2	Fujifilm 500/600	10	FAL23607
MDD	WB-OL2-NC2	Olympus	10	FAL23673
MDD	WB-PX4-WC2	Pentax	10	FAL23470

## HybriX EcoFlushCaps with CO2



- Reduced risk of cross-contamination by eliminating the need to clean and reprocess reusable sterile water bottles
- The sterile system can be used for 24 hours
- For connection to all common flushing pump units
- With CO2 hose system
- With hose set and screw cap for all common sterile water bottles
- Flushing including check valve

Supplier	Product Code	Compatibility	QOM	FAL Code
MDD	WB-HYB-FJ4	Fujifilm 500/600	10	FAL24554
MDD	WB-HYB-FJ5	Fujifilm 700	10	FAL24583
MDD	WB-HYB-OL4	Olympus	10	FAL24597
MDD	WB-HYB-PX5	Pentax	10	FAL24609



- Clearly labelled connection points allow the end user to quickly set up the tubing and prepare the endoscope and theatre

Supplier	Product Code	Description	Compatibility	QOM	FAL Code
MDD	IT-ETS-210	Irrigation Tube With Luer-Lock Connector Including Check-Valve	All Common Sterile Water Bottles	15	
MDD	IT-OUS-210	Irrigation Tube	Olympus, Erbe & Medivator Pumps	15	FAL23584
<b>CO2 Device Adapters</b>					
MDD	ADA-CO2-FJ	CO2 Adapter	Fujifilm	100	FAL24559
MDD	ADA-CO2-OUS	CO2 Adapter	Olympus	100	FAL24610



- We offer a full complement of CO2 insufflation tubing compatible with most main line pump manufacturers
- Clearly labelled connection points allow the end user to quickly set up the tubing and prepare the endoscope and theatre

Supplier	Product Code	Description	Compatibility	QOM	FAL Code
Vytil	VTL/CN04	Connector Tube - Re-Usable	Olympus	20	
Vytil	VTL/TB0403EXT	Inflation Tube With CO2	Pentax	10	
Vytil	VTL/TB0406	Inflation Tube	Pentax	10	
Vytil	VTL/TB1101	Inflation Tube With CO2	Fujifilm 500/600	10	
Vytil	VTL/TB1102EXT	Inflation Tube With CO2	Fujifilm 500/600	10	
Vytil	VTL/TB1201	Inflation Tube With CO2	Fujifilm 700	10	
Vytil	VTL/TB1202EXT	Inflation Tube With CO2	Fujifilm 700	10	
MDD	CO2-TU-210	CO2 Connection Tube 120cm		10	FAL23738
MDD	OLY-CT-001	Single Use Cleaning Tube	Olympus	50	FAL24566



Olympus Adapter



Pentax Adapter



Fujifilm Adapter

- Single use Jet adapter option for endoscopes offering an auxiliary water channel
- Easy to identify design and fixture allows for ease of set up and removes cross contamination risks associated with reusable devices

Supplier	Product Code	Compatibility	QOM	FAL Code
MDD	ADA-02-JC-FJ	Fujifilm 700	100	FAL24605
MDD	ADA-02-JC-OUS	Olympus Plastic	100	FAL24578
MDD	ADA-02-JC-PE	Pentax Plastic	100	FAL24590
MDD	ADA-JC-FJ	Fujifilm 700	100	FAL23663
MDD	ADA-JC-OUS	Olympus	100	FAL23733
MDD	ADA-JC-PE	Pentax	100	



- Lubricating water soluble jelly, for endoscopy, proctology and more in general

Supplier	Product Code	Description	Volume (g)	QOM	FAL Code
Vytil	LUB0005	Sterile Lubricating Jelly - Foil Packet	5	150	
Vytil	LUB0082	Sterile Lubricating Jelly - Flip Top Tube	82	1	

# ENDO THERAPY







Lugol life 2%



Acetic Acid



Carmine Indigo

### Lugol life 2%

- In digestive endoscopy: demonstration of epidermoid lesions of the oesophagus
- In colposcopy: brushing the cervix - Schiller test


### Carmine Indigo

- Enhancement of mucosal relief anomalies by spraying on the colonic wall
- Recommended concentration for use from 0.2% to 0.4%
- Not sterile, not injectable

### Acetic Acid


- Mucosa abnormalities enhancement agent
- Non-sterile solution of acetic acid 3% or 5%
- Use within 30 days after opening

Supplier	Product Code	Description	Concentration	QOM	FAL Code
Vytil	849	Lugol Life Glass Ampoule of 10 ml	2%	10	FAL23556
Vytil	J2F/810160/4	Indigo Carmine Glass Bottle of 10 ml	1%	10	
Vytil	ACE3100	Acetic acid 100 ml bottle	3%	20	FAL23835
Vytil	ACE5100	Acetic acid 100 ml bottle	5%	20	ACE5100




- Used in conjunction with an endoscope to spray the gastrointestinal mucosa
- Metal tip spray catheter for accurate 360° forward spray action for even and precise dying of the mucosa
- Ideal for use with diluted Indigo carmine for chromoendoscopy surveillance
- Available in a 220cm length for both gastroscopy and colonoscopy

Supplier	Product Code	Description	Catheter diameter (mm)	Usable length (cm)	Operating channel diameter (mm)	QOM	FAL Code
Vytil	NLS/SP-18-220	Spray Catheter	1.8	220	> 2.0	10	FAL23605
Vytil	NLS/SP-23-230	Spray Catheter	2.3	230	> 2.8	20	FAL24600
Vytil	NLS/SP-18-120	Spray Catheter	1.8	120	> 2.0	10	FAL23885
Pentax	ST-23-230-A	Spray Catheter - Straight Spray	2.3	230	> 2.8	10	
Pentax	ST-23-230-B	Spray Catheter - Water Fog Spray	2.3	230	> 2.8	10	



- Available for oesophagus, cardia, pylorus, colon and biliary duct.
- All balloons are TTS compatible except achalasia balloons
- One size fits all
- 2 radiopaque marks
- Compatible with guidewires of 0.035 "maximum

Supplier	Product Code	Balloon		Catheter		Pressure (ATM) Recommended / Burst	Comments	QOM	FAL Code
		Diameter (mm)	Length (mm)	Diameter (mm/fr)	Length (cm)				
<b>Dilation balloon single diameter</b>									
Vytil	CJ-JHY-BD-06-40-250	6	40	2.4 / 7	250	8 - 10			FAL24573
Vytil	CJ-JHY-BD-08-40-250	8	40	2.4 / 7	250	8 - 10			FAL24593
Vytil	CJ-JHY-BD-10-40-250	10	40	2.4 / 7	250	8 - 10	Operating channel diameter: ≥ 2.8 mm	1	FAL24577
Vytil	CJ-JHY-BD-12-60-250	12	60	2.4 / 7	250	5 - 6			FAL24568
Vytil	CJ-JHY-BD-14-60-250	14	60	2.4 / 7	250	5 - 6	Guide wire: 0.035"		FAL24557
Vytil	CJ-JHY-BD-16-60-250	16	60	2.4 / 7	250	5 - 6		FAL24553	
Vytil	CJ-JHY-BD-16-80-250	16	80	2.4 / 7	250	5 - 6		FAL24550	
Vytil	CJ-JHY-BD-18-80-250	18	80	2.4 / 7	250	4 - 5			FAL24580
Vytil	CJ-JHY-BD-20-80-250	20	80	2.4 / 7	250	4 - 5			FAL24552
<b>Achalasia balloon single diameter</b>									
Vytil	CJ-JHY-BD-30-80-90	30	80	6.7 / 20	90	2 - 3	OTW procedure guide wire: 0.035"	1	FAL24596
Vytil	CJ-JHY-BD-35-80-90	35	80	6.7 / 20	90	2 - 3			FAL24562
Vytil	CJ-JHY-BD-40-80-90	40	80	6.7 / 20	90	2 - 3			FAL24565



- Soft and tapered tip for smooth advancement through the stricture
- 3 Size Diameter Balloons
- Compatible with inflation gun type device & fluxy gun type

Supplier	Product Code	Balloon		Catheter		Pressure (ATM) Recommended / Burst	Comments	QOM	FAL Code
		Diameter (mm)	Length (mm)	Diameter (mm/fr)	Length (cm)				
<b>Compliant dilation balloon 3 diameters pre-mounted on 0.035" guide wire</b> Guide length: 300cm / hydrophilic length: 5.5cm									
Vytil	NLS/BD3-8-55-240-GW	8 / 9 / 10	55	2.4	240	3 / 4.5 / 6	Including guide wire : 0.035"	1	FAL24572
Vytil	NLS/BD3-10-55-240-GW	10 / 11 / 12	55						FAL24571
Vytil	NLS/BD3-12-55-240-GW	12 / 13.5 / 15	55						FAL24579
Vytil	NLS/BD3-15-55-240-GW	15 / 16.5 / 18	55						FAL24589
Vytil	NLS/BD3-18-55-240-GW	18 / 19 / 20	55						FAL24560
Vytil	NLS/BD3-12-80-240-GW	12 / 13.5 / 15	80						FAL24591
Vytil	NLS/BD3-15-80-240-GW	15 / 16.5 / 18	80						FAL24581
Vytil	NLS/BD3-18-80-240-GW	18 / 19 / 20	80						FAL24603



- Endoscopic dilation bougie, with radio-opaque markers and graduation
- Wide range of diameter to meet different needs.

Supplier	Product Code	Description	Diameter (mm)	QOM	FAL Code
Vytil	CJ-BOUGIES5	Disposable Bougie 85 cm length	5	1	
Vytil	CJ-BOUGIES6		6		
Vytil	CJ-BOUGIES7		7		
Vytil	CJ-BOUGIES8		8		
Vytil	CJ-BOUGIES9		9		
Vytil	CJ-BOUGIES10		10		
Vytil	CJ-BOUGIES11		11		
Vytil	CJ-BOUGIES12		12		
Vytil	CJ-BOUGIES13		13		
Vytil	CJ-BOUGIES14		14		
Vytil	CJ-BOUGIES15		15		
Vytil	CJ-BOUGIES16		16		
Vytil	CJ-BOUGIES17		17		
Vytil	CJ-BOUGIES18		18		
Vytil	CJ-BOUGIES19		19		
Vytil	CJ-BOUGIES20		20		



- Progressive Balloon Dilatation Catheter provides consistent performance for balloon endoscopy with optimal control, efficiency and performance.
- Amara Wire-guided Progressive Balloon Dilatation Catheters are also indicated for use in removal of difficult biliary stones (DASE)
- Made with a new balloon folding process for lower profile leading to consistent ease of passage through the scope
- Embedded Radiopaque markers to facilitate balloon placement using fluoroscopy
- Supplied with a 0.035" Guidewire, clamp and stopcock to facilitate procedure time

Supplier	Product Code	Balloon		Catheter		Inflation Pressure (ATM)	Comments	QOM	FAL Code
		Inflated Diameter (mm)	Length (mm)	Diameter (mm/fr)	Length (cm)				
Envaste	GDB-0110655S	6 / 7 / 8	5.5	7.4 / 2.45	180	3 / 6 / 10	Including guide wire : 0.035"	1	FAL23653
Envaste	GDB-0110855S	8 / 9 / 10				3 / 5.5 / 9			FAL23498
Envaste	GDB-0111055S	10 / 11 / 12				3 / 5 / 8			FAL23602
Envaste	GDB-0111255S	12 / 13.5 / 15				3 / 4.5 / 8			FAL23477
Envaste	GDB-0111555S	15 / 16.5 / 18				3 / 4.5 / 7			FAL23644
Envaste	GDB-0111855S	18 / 19 / 20			3 / 4.5 / 6	FAL23848			
Envaste	GDB-0120655S	6 / 7 / 8			3 / 6 / 10	FAL23758			
Envaste	GDB-0120855S	8 / 9 / 10			3 / 5.5 / 9	FAL23476			
Envaste	GDB-0121055S	10 / 11 / 12			3 / 5 / 8	FAL23691			
Envaste	GDB-0121255S	12 / 13.5 / 15			3 / 4.5 / 8	FAL23874			
Envaste	GDB-0121555S	15 / 16.5 / 18	3 / 4.5 / 7	FAL23834					
Envaste	GDB-0121855S	18 / 19 / 20	3 / 4.5 / 6	FAL23609					




- Used to supply pressure through the catheter in order to inflate or deflate the balloon
- Compatible with any type of Luer Lock connection

Supplier	Product Code	Description	Maximal pressure (ATM / PSI)	Maximum volume (ml)	QOM	FAL Code
Vytil	NLS/ID20G	Gun Type	30 / 440	20	20	FAL23480
Vytil	NLS/ID20A	Fluxy	30 / 440	30	20	



- Slit and Hole - Fluid in the hood is discharged by capillarity and pressure difference

Supplier	Product Code	Hood Protrusion Length (mm)	Endoscope Diameter Compatibility (mm)	QOM	FAL Code
Pentax	SHF-010	3.5 ~ 7.0	9.0 ~ 9.2	5	
Pentax	SHF-025		9.8 ~ 9.9	5	
Pentax	SHF-030		9.9 ~ 10.5	5	
Pentax	SHF-035		10.8	5	
Pentax	SHF-060		11.7 ~ 12.2	5	
Pentax	SHF-050		12.8	5	
Pentax	SHF-055		13.2	5	



### Superior Navigation

- Inflatable two-balloon design stabilises the endoscope and the nearby anatomy
- Intraoperative conduit enables safe and efficient access, and repeated removal of tools and tissue, without losing position

### Enhanced Visualisation

- Aft balloon supplies gentle endoluminal positioning, helping to reveal tissue planes and expose hidden polyps within folds and flexures.
- EZ Glide hydrophilic coating provides a clearer, more visible therapeutic zone without messy lubricants

### Stability & Control

- Colon shortening allows for safe and gentle straightening of tortuous anatomy, achieving faster access and fewer aborted procedures
- Extendable fore balloon creates a stable therapeutic zone, flattening the surgical site for submucosal injection and precise tissue dissection
- Integrated suture system supplies controlled, dynamic tissue manipulation and tension during procedures

Supplier	Product Code	Length (cm)	QOM	FAL Code
Lumendi	D-1001-5	130	5	
Lumendi	D-1002-5	103	5	
Lumendi	D-1001	130	1	
Lumendi	D-1002	103	1	

- Disk in the knife helps to hook up the tissue for improved safety and ease of cutting
- Metal plate that closes the knife for clear marking and better hemostasis
- Water jet function that enables additional injection without changing the device

Supplier	Product Code	Knife		Max. Insertion Portion Width (mm)	Catheter		QOM	FAL Code
		Knife Length Extended (mm)	Knife Length Retracted (mm)		Diameter (mm)	Length (cm)		
Pentax	DN-D2718B	2	0.5	2.7	2.8	180	1	
Pentax	DN-D2722B	2	0.5	2.7	2.8	220	1	




- Insulated tip and single-sided knife designed for a safer dissection
- Rotatable tip with optical marker supports a precisely directed dissection
- Extended Knife length designed for a faster procedure

Supplier	Product Code	Max. Insertion Portion Width (mm)	Knife Length (mm)	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
Pentax	DP-D2518	2.5	5	2.8	180	1	
Pentax	DP-D2622	2.6	2.5	2.8	220	1	



- Bipolar HemoStat WideCup reduces the impact of deeper tissue thermal degeneration
- Supports prevention of re-bleeding
- Based on its flexibility, it can be used for various applications such as EMR/ESD, gastrointestinal bleeding and papilla resection

Supplier	Product Code	WideCup Opening (mm)	Max. Insertion Portion Width (mm)	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
Pentax	HS-D2618	5.8	2.6	2.8	180	5	
Pentax	HS-D2622	5.8	2.6	2.8	220	5	




• DiLumen +DRTM Dynamic Retraction System consists of a controllable sliding balloon assembly with two suture loops of 2 cm and 5 cm diameters. Depending on polyp position, either loop can be fixed to mucosal tissue via a common hemostasis clip and pushed or pulled for improved tissue exposure

### DiLumen +DR provides:

- Effective retraction of tissue
- Continuous tissue tension
- Visualization of submucosal tissue plane
- Precise dissection

Supplier	Product Code	Length (cm)	QOM	FAL Code
Lumendi	D-1001-5	130	5	
Lumendi	D-1002-5	103	5	
Lumendi	D-1001	130	1	
Lumendi	D-1002	103	1	



- Profiled tip for a precise insertion
- Rotatable, visual control on the handle
- Triple lumen
- Compatible short and long guidewire

Supplier	Product Code	Description	Extremity length (cm)	Cutting wire length (mm)	Catheter length (cm)	QOM	FAL Code
Vytil	NLS/STU-3RX-05-25	Sphincterotome for short and long guidewire	5	25	200	5	FAL23744
Vytil	NLS/STU-3RX-05-30	Sphincterotome for short and long guidewire	5	30	200	5	FAL23905
Vytil	NLS/ST-3RX-05-25	Sphincterotome PTFE covered for short and long guidewire	5	25	200	5	FAL23667
Vytil	NLS/ST-3RX-05-30	Sphincterotome PTFE covered for short and long guidewire	5	30	200	5	FAL23677



### Stone Extraction Basket - Rotatable and wire-guided

- Nitinol wires with shape memory (4 single-strand)
- Rotary spiral shape
- Compatible with short or long guide wire 0.035" maximum


Supplier	Product Code	Basket diameter (mm)	Catheter Length (cm)	Diameter (mm)	Minimal operating channel diameter (mm)	QOM	FAL Code
Vytil	NLS/BKS-20-24-195G	20	195	3.1	3.7	5	FAL23620
Vytil	NLS/BKS-30-24-195G	30	195	3.1	3.7	5	FAL23725



- Stone retrieval balloon, inflation with air only
- Compatible long and short guide wires
- 3 size balloon, triple lumen
- Radiopaque markers

Supplier	Product Code	Balloon Diameter (mm)	Length Catheter (cm)	Diameter (mm)	Compatible Guide Wire	QOM	FAL Code
Vytil	NLS/BL-15-23-200	10	200	2.3	0.035"	5	FAL23786
		12					
		15					
Vytil	NLS/BL-L-23-200	13	200	2.3	0.035"	5	FAL23897
		15					
		18					





- Storage of accessories in the side pockets.
- Ideal for storing devices and supplies that can inadvertently fall on the floor, get misplaced, tangled, or damaged during procedures.


Supplier	Product Code	Description	Dimensions (cm)	QOM	FAL Code
Vytil	NLS/SF3-M	Sterile surgical drape with 3 pockets (Kangoo)	75 x 130	100	FAL24594



- Innovative handle enabling a 360° rotation with one hand
- Repositionable
- Optimal flexibility of the distal part
- Short tip after release: 7 mm (up to 13 mm max. open) to 9 mm (16 mm max. open)
- MRI compatible

Supplier	Product Code	Maximal Opening (mm)	Jaws Length (mm)	Jaws Angle	Diameter (mm)	Length (cm)	Handle Colour	QOM	FAL Code
Vytil	CJ-HC09-25-230	9	5	135°	2.5	230	■	10	FAL24588
Vytil	CJ-HC11-25-160	11	6	135°	2.5	160	■		
Vytil	CJ-HC11-25-230	11	6	135°	2.5	230	■		
Vytil	CJ-HC11-25-260	11	6	135°	2.5	260	■		
Vytil	CJ-HC13-25-230	13	7	135°	2.5	230	■		FAL23515
Vytil	CJ-HC1690-25-230	16	9	90°	2.5	230	■		FAL23559

# NovaClip R4



- Innovative handle enabling a 360° rotation with one hand
- Repositionable
- Short tip after release: 7 mm (12 clip) and 8 mm (15 clip)
- Optimal flexibility of the distal part
- MRI compatible

Supplier	Product Code	Maximal Opening (mm)	Jaws Angle	Diameter (mm)	Length (cm)	Handle Colour	QOM	FAL Code
Vytal	VTL/HR-230-12-135-Y	12	135°	2.6	230	■	10	FAL23635
Vytal	VTL/HR-230-15-135-Y	15	135°	2.5	230	■	10	FAL23624

# HemoClips




- Synchronous rotation supports accurate positioning
- Repositionable clips
- The handle's security lock system helps prevent unexpected pain for the patient
- MR conditional

Supplier	Product Code	Maximal Opening (mm)	Catheter		QOM	FAL Code
			Diameter (mm)	Length (cm)		
Pentax	CD-11-230-O	11	2.8	230	5	
Pentax	CD-11-230-P	11	2.8	230	5	
Pentax	CD-13-230-P	13	2.8	230	5	
Pentax	CD-16-230-P	16	2.8	230	5	

# Novaject EMR Needles

## Injection Needles



Triple bevel needle for an easy penetration

- Ergonomic handle to control needle movements
- Transparent catheter for liquid flow visualisation

Supplier	Product Code	Needle		Catheter		QOM	FAL Code
		Diameter (G / mm)	Length (mm)	Diameter (mm)	Length (cm)		
Vytil	CJ-A3-23-230-21	21 / 0.8	3	2.3	230	10	FAL24588
Vytil	CJ-A3-23-230-23	23 / 0.6	3	2.3	230		FAL24563
Vytil	CJ-A3-23-230-25	25 / 0.5	3	2.3	230		FAL24576
Vytil	CJ-A-18-165-25	25 / 0.5	4	1.8	165		FAL23879
Vytil	CJ-A-23-160-19	19 / 1.1	4	2.3	160		FAL23510
Vytil	CJ-A-23-230-19	19 / 1.1	4	2.3	230		FAL23687
Vytil	CJ-A-23-230-21	21 / 0.8	4	2.3	230		FAL23507
Vytil	CJ-A-23-230-23	23 / 0.6	4	2.3	230		FAL23864
Vytil	CJ-A-23-230-25	25 / 0.5	4	2.3	230		FAL23642
Vytil	CJ-A5-23-230-21	21 / 0.8	5	2.3	230		FAL23850
Vytil	CJ-A5-23-230-23	23 / 0.6	5	2.3	230		FAL23488
Vytil	CJ-A5-23-230-25	25 / 0.5	5	2.3	230		FAL23729

# Injection Needles



- Ergonomic handle to control needle movements
- Transparent catheter for liquid flow visualisation

Supplier	Product Code	Needle		Catheter		QOM	FAL Code
		Diameter (G / mm)	Length (mm)	Diameter (mm)	Length (cm)		
MDD	IN-23-19-4-230	19 / 1.0	4	2.4	230	10	FAL23697
MDD	IN-23-19-5-230	19 / 1.0	5	2.4	230		
MDD	IN-23-19-6-230	19 / 1.0	6	2.4	230		FAL23656
MDD	IN-23-22-4-160	22 / 0.7	4	2.4	160		
MDD	IN-23-22-4-230	22 / 0.7	4	2.4	230		FAL23521
MDD	IN-23-22-5-160	22 / 0.7	5	2.4	160		
MDD	IN-23-22-5-230	22 / 0.7	5	2.4	230		
MDD	IN-23-22-6-230	22 / 0.7	6	2.4	230		FAL23517
MDD	IN-23-23-4-160	23 / 0.6	4	2.4	160		
MDD	IN-23-23-4-230	23 / 0.6	4	2.4	230		FAL23831
MDD	IN-23-23-5-160	23 / 0.6	5	2.4	160		
MDD	IN-23-23-5-230	23 / 0.6	5	2.4	230		FAL23679
MDD	IN-23-23-6-230	23 / 0.6	6	2.4	230		FAL23574
MDD	IN-23-25-4-230	25 / 0.5	4	2.4	230		FAL23898
MDD	IN-23-25-5-230	25 / 0.5	5	2.4	230		FAL23630
MDD	IN-23-25-6-230	25 / 0.5	6	2.4	230	FAL23572	
Pentax	IN-18-120-2205-A	22 / 0.5	120	1.8	180		
Pentax	IN-23-180-2204-A	22 / 0.4	180	2.4	180		
Pentax	IN-23-230-1904-A	19 / 0.4	230	2.4	230		
Pentax	IN-23-230-2205-A	22 / 0.5	230	2.4	230		
Pentax	IN-23-230-2504-A	25 / 0.4	230	2.4	230		
Pentax	IN-23-230-2506-A	25 / 0.6	230	2.4	230		

# Cold Snares


## Snares



- Cold snare for small polyps resection
- Thin and stiff braided wire
- Rotatable or not


Supplier	Product Code	Description	Snare Width (mm)	Snare Diameter (mm)	Catheter		QOM	FAL Code
					Diameter (mm)	Length (cm)		
Vytil	CJ-CS-23-230-D10	Cold Snare	10	0.28	2.3	230	10	FAL24606
Vytil	CJ-CSR-23-230-D10	Rotatable Cold Snare	15	0.28	2.3	230		FAL24555
MDD	MST-CR-24-02-10-230	Rotatable Cold Snare	10	0.25	2.4	230		FAL23581
MDD	MST-CR-24-03-10-230	Rotatable Cold Snare	10	0.30	2.4	230		FAL23492
MDD	MST-CS-24-02-10-230	Cold Snare	10	0.25	2.4	230		FAL23698
MDD	MST-CS-24-03-10-230	Cold Snare	10	0.30	2.4	230		FAL23474
Pentax	CS-23-230-10-C3	Cold Snare	10		2.3	230		
Pentax	CS-23-230-15-C3	Cold Snare	15		2.3	230		

# Rotatable Snares



Supplier	Product Code	Snare Width (mm)	Catheter		QOM	FAL Code
			Diameter (mm)	Length (cm)		
Vytil	CJ-ADR2-23-230-O10	10	2.3	230	10	FAL23681
Vytil	CJ-ADR2-23-230-O15	15	2.3	230		FAL23631
Vytil	CJ-ADR2-23-230-O20	20	2.3	230		FAL23672
Vytil	CJ-ADR2-23-230-O25	25	2.3	230		FAL23901
Vytil	CJ-ADR2-23-230-O30	30	2.3	230		FAL23803
Vytil	CJ-ADR2-23-300-O30	30	2.3	300		FAL23828
Vytil	CJ-ADR-23-230-O25	25	2.3	230		FAL23588
Vytil	VTL/SN-24-240-O15R	15	2.4	240		FAL23717
Vytil	VTL/SN-24-240-O20R	20	2.4	240		FAL23603
Vytil	VTL/SN-24-240-O25R	25	2.4	240		FAL23707
Vytil	VTL/SN-24-240-O30R	30	2.4	240	FAL23871	
MDD	OVA-DR-24-10-230	10	2.4	230	10	FAL24587
MDD	OVA-DR-24-15-160	15	2.4	160		FAL23668
MDD	OVA-DR-24-15-230	15	2.4	230		FAL23801
MDD	OVA-DR-24-20-160	20	2.4	160		FAL23552
MDD	OVA-DR-24-20-230	20	2.4	230		FAL23542
MDD	OVA-DR-24-25-160	25	2.4	160		FAL23875
MDD	OVA-DR-24-25-230	25	2.4	230		FAL23773
MDD	OVA-DR-24-30-160	30	2.4	160		FAL23482
MDD	OVA-DR-24-30-230	30	2.4	230		FAL23479
MDD	OVA-DR-24-50-230	50	2.4	230		FAL23689
Pentax	SD-23-230-20-A1R	20	2.3	230		
Pentax	SD-23-230-30-A1R	30	2.3	230		

# Monofilament Snares



CJ-ADR2-23-230-O25B




CJ-ADR-23-230-MFO-O20

- Stiff single element monofilament wire provides maximum snare rigidity for challenging EMR
- Highly effective on scarred or fibrotic tissue
- Unique Distal Traction Pin increases snare durability along with top anchoring
- Clear white ergonomic handle for clear marking
- Diameter wire 0.3 mm

Supplier	Product Code	Snare Width (mm)	Description	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
Vytil	CJ-ADR-23-230-MFO10	10	Hot & Cold Snare	2.3	230	10	
Vytil	CJ-ADR-23-230-MFO15	15	Hot & Cold Snare	2.3	230		
Vytil	CJ-ADR-23-230-MFO20	20	Standard EMR Snare	2.3	230		
Vytil	CJ-ADR2-23-230-O25B	15 / 25	Two in One	2.3	230		FAL23626

# Multi-Stage Snares



360°

- Monofilament snares: 0.3 mm single wire
- Braided snares: 0.35 or 0.36 mm braided stiff or rotatable wire

Supplier	Product Code	Snare Width (mm)	Rotatable	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
MDD	OVA-HYB-24-10-230	10	Yes	2.4	230	10	FAL23845
MDD	OVA-HYB-24-15-230	15	Yes	2.4	230		FAL23746

Supplier	Product Code	Description	Snare Width (mm)	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
Vytil	CJ-ADR-23-230-O10	Oval	10 x 15	2.3	230	10	
Vytil	CJ-ADR-23-230-O15		15 x 20	2.3	230		
Vytil	CJ-ADR-23-230-O20		20 x 25	2.3	230		
Vytil	CJ-ADR-23-230-O30		30 x 45	2.3	230		FAL23655
Vytil	CJ-ADR-23-230-O35		35 x 50	2.3	230		FAL23616
Vytil	CJ-ADR-23-230-O40		40 x 60	2.3	230		FAL23829
Vytil	CJ-ADR-23-230-O45		45 x 65	2.3	230		FAL23516
Vytil	CJ-ADR-23-230-O50		50 x 70	2.3	230		FAL23890
Vytil	VTL/SN-24-240-O10		10	2.4	240		FAL23730
Vytil	VTL/SN-24-240-O15		15	2.4	240		FAL23501
Vytil	VTL/SN-24-240-O20		20	2.4	240		FAL23595
Vytil	VTL/SN-24-240-O25		25	2.4	240		FAL23596
Vytil	VTL/SN-24-240-O30		30	2.4	240		FAL23699
MDD	OVA-DS-23-10-230		10	2.3	230		FAL23823
MDD	OVA-DS-23-15-160		15	2.3	160		FAL23540
MDD	OVA-DS-23-15-230		15	2.3	230		FAL23641
MDD	OVA-DS-23-20-160		20	2.3	160		FAL23768
MDD	OVA-DS-23-20-230		20	2.3	230		FAL23892
MDD	OVA-DS-23-25-160		25	2.3	160		FAL23524
MDD	OVA-DS-23-25-230	25	2.3	230	FAL23783		
MDD	OVA-DS-23-30-160	30	2.3	160	FAL23680		
MDD	OVA-DS-23-30-230	30	2.3	230	FAL23841		
MDD	OVA-DS-23-35-230	35	2.3	230	FAL24586		
MDD	OVA-DS-23-50-230	35	2.3	230	FAL23742		
Pentax	SD-23-230-15-A1	15	2.3	230			
Pentax	SD-23-230-25-A1	25	2.3	230			
Pentax	SD-23-230-35-A1	35	2.3	230			
Pentax	SD-23-230-50-A1	50	2.3	230			



Hexagonal

- The memory points and braided wire features of our hexagonal snare provides stiffness and reduction in shape deformity

Supplier	Product Code	Description	Snare Width (mm)	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
MDD	HEX-DS-23-10-230	Hexagonal	10	2.4	230	10	FAL23911
MDD	HEX-DS-23-15-160		15	2.4	160		FAL23577
MDD	HEX-DS-23-15-230		15	2.4	230		FAL23852
MDD	HEX-DS-23-20-160		20	2.4	160		FAL23772
MDD	HEX-DS-23-20-230		20	2.4	230		FAL23551
MDD	HEX-DS-23-25-160		25	2.4	160		FAL23855
MDD	HEX-DS-23-25-230		25	2.4	230		FAL23705
MDD	HEX-DS-23-30-160		30	2.4	160		FAL23856
MDD	HEX-DS-23-30-230		30	2.4	230		FAL23870
MDD	HEX-DS-23-35-230		35	2.4	230		FAL24601
Vytil	CJ-ADR-23-230-H15		15	2.3	230		FAL23825
Vytil	CJ-ADR-23-230-H25		25	2.3	230		FAL23757
Vytil	CJ-ADR-23-230-H30		30	2.3	230		FAL23485
Vytil	CJ-ADR-23-230-H35		35	2.3	230		FAL23849
Pentax	SD-18-160-20-C1	20	1.8	160			



Crescent

- The memory points and braided wire features of our hexagonal snare provides stiffness and reduction in shape deformity

Supplier	Product Code	Description	Snare Width (mm)	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
MDD	CRE-DS-23-10-230	Crescent Snares	10	2.3	230	10	FAL23794
MDD	CRE-DS-23-15-160		15	2.3	160		FAL23632
MDD	CRE-DS-23-15-230		15	2.3	230		FAL23503
MDD	CRE-DS-23-20-160		20	2.3	160		FAL23756
MDD	CRE-DS-23-20-230		20	2.3	230		FAL23600
MDD	CRE-DS-23-25-160		25	2.3	160		FAL23878
MDD	CRE-DS-23-25-230		25	2.3	230		FAL23508
MDD	CRE-DS-23-30-160		30	2.3	160		FAL23816
MDD	CRE-DS-23-30-230		30	2.3	230		FAL23882
MDD	CRE-DS-23-35-230		35	2.3	230		FAL24570
Vytil	CJ-ADR-23-230-C15		15	2.3	230		FAL23868
Vytil	CJ-ADR-23-230-C25		25	2.3	230		FAL23601




- Round snare loop design aides precise approximation of lesion by mirroring typical small (<15mm) lesion pathology
- Round loop shape shifts emphasis to lateral wire closure, minimising distal loop slipping, increasing control and effective resection

Supplier	Product Code	Description	Snare Width (mm)	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
Vytil	CJ-ADR-23-230-HCO10	Oval Snares	10	2.3	230	10	FAL23652
Vytil	CJ-ADR-23-230-HCO15		15	2.3	230		FAL23541
Vytil	NLS/PS-23-230-O10		10	2.4	230		FAL23666
Vytil	NLS/PS-23-230-O15	Hybrid Hot/Cold Snares	15	2.4	230	10	FAL23724
Vytil	VTL/SN-24-240-HCO15	Rotatable Hybrid Hot/Cold Snares	15	2.4	240	10	FAL23716
Vytil	VTL/SN-24-240-HCO10R		10	2.4	240		FAL23569
Vytil	VTL/SN-24-240-HCO15R		15	2.4	240		FAL23649



Spiral snare

- Spiral snare has additional spiral wire around the snare to prevent slipping on the mucosa
- Spiral snare has an ultra stiff wire
- Also a standard hot rotatable snare is available




- Multi-functioning design helps address a variety of clinical needs
- Three-shapes in one device can help towards a wide range of clinical requirements when encountered
- 3 distinct shapes and sizes should eliminate the need to exchange snares
- Ideal for use when facing varying clinical challenges
- The stiffer braided wire minimises distal deflection

Supplier	Product Code	Snare Width (mm)	Catheter		QOM	FAL Code
			Diameter (mm)	Length (cm)		
Vytil	CJ-ADR-23-230-HCO06	6	2.3	230	10	
Vytil	VTL/SN-24-240-HCO10	10	2.3	230		FAL23678
Vytil	CJ-ADR2-23-230-NG20	20 x 25	2.3	230		FAL23619
Vytil	CJ-ADR2-23-230-NG25	25 x 35	2.3	230		FAL23782
Vytil	CJ-ADR2-23-230-NG30	30 x 45	2.3	230		FAL23889


Supplier	Product Code	Description	Snare Width (mm)	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
MDD	MOE-24-6-12-30-230	Standard Multi-Stage Snare	30 / 10 / 6	2.4	230	10	FAL24575
MDD	MOS-24-6-12-30-230	Rotatable Multi-Stage Snare	30 / 10 / 6	2.4	230		FAL24595
MDD	MOT-24-6-12-30-230	Teflon Coated Multi-Stage Snare	30 / 10 / 6	2.4	230		FAL24608





- Enables hybrid techniques EMR/ESD especially in the colon
- Combines an adjustable knife on the tip designed for precise marking and incision, and a snare for resection
- The all-in-one accessory for supporting a complete lesion resection

Supplier	Product Code	Snare Width (mm)	Knife Length (mm)	Adjustable Knife Length (mm)	Catheter		QOM	FAL Code
					Diameter (mm)	Length (cm)		
Pentax	DO-D2618	22	6	2 - 6	2.8	180	1	
Pentax	DO-D2622	22	6	2 - 6	2.8	220		



- Rotatable retrieval net
- Foreign body and polyp retrieval

Supplier	Product Code	Snare width (mm)	Description	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
Vytil	VTL/RRN-18-160-15	15	Standard Oval	1.8	160	10	
Vytil	VTL/RRN-18-160-20	20	Standard Oval	1.8	160		
Vytil	VTL/RRN-23-230-30	30	Standard Oval	2.3	230		
Vytil	VTL/RRN-23-230-40	40	Standard Oval	2.3	230		
Vytil	VTL/RRN-23-300-40	40	Standard Oval	2.3	300		
MDD	FBN-O35-25-230-SU	35 x 70	Octagonal	2.5	230		FAL23836
MDD	FBN-OR35-25-230-SU	35 x 70	Octagonal	2.5	230		FAL23575
MDD	FBN-S30-25-230-SU	20 x 30	Standard Oval	2.5	230		FAL23796
MDD	FBN-S50-25-230-SU	25 x 50	Standard Oval	2.5	230		FAL23627
MDD	FBN-S60-25-230-SU	30 x 60	Standard Oval	2.5	230		FAL23910
MDD	FBN-S75-25-230-SU	40 x 70	Standard Oval	2.5	230	FAL23472	
Pentax	FG-25-230-C3	30 x 60	Standard Oval	2.5	230		
Pentax	FG-25-230-C4	40 x 75	Standard Oval	2.5	230		
Pentax	FG-25-230-C5	40 x 65	Standard Oval	2.5	230		



PT-01-05



VTL/PT-02

### Two Drawer Polyp Trap

- The soft catch tray system permits effective retrieval of pathology without compromising suction and loss of procedure time
- Integrity of pathology maintained with soft catch drawer design

### Five Drawer Polyp Trap

- 5 clearly marked chambers to match polyp excision location
- Transparent chambers enables identification of polyp

Supplier	Product Code	Description	QOM	FAL Code
Vytil	VTL/PT-01	One Chamber Polyp Trap	25	FAL23721
Vytil	VTL/PT-02	Two Chamber Polyp Trap	10	FAL24652
Vytil	VTL/PT-04	Four Chamber Polyp Trap	25	
Vytil	VTL/PT-05		25	FAL23660
MDD	PT-01-05	Five Chamber Polyp Trap	10	FAL23720
Pentax	JHY-PT-05		10	



Alligator

Easy Grasper

Pelican

Prong Grasper

- Alligator: excellent grip force avoiding the deterioration of the sample
- Easy Grasper: for an excellent grasping
- Pelican: excellent grip force avoiding loss of the sample
- Prong Grasper: Atraumatic prongs securely capture pathology without damaging the sample integrity

Supplier	Product Code	Description	Catheter (mm)	Length (cm)	QOM	FAL Code
Vytil	CJ-2418GH	Alligator	2.3	230	5	FAL23853
Vytil	CJ-2423GK	Pelican	2.3	230	5	FAL23659
Vytil	CJ-PPDR-23-180	Rat Tooth	2.3	180	5	FAL23812
Vytil	CJ-PPF-23-230	4 Prong	2.3	230	5	FAL23547
MDD	GFR-ALI-16-61-SU	Rat Tooth Alligator	1.6	61	10	FAL24574
MDD	GFR-ALI-18-120-SU	Rat Tooth Alligator	1.8	120	10	FAL23784
MDD	GFR-ALI-18-40-SU	Rat Tooth Alligator	1.8	40	10	FAL24598
MDD	GFR-ALI-23-160-SU	Rat Tooth Alligator	2.3	160	10	FAL23832
MDD	GFR-ALI-23-230-SU	Rat Tooth Alligator	2.3	230	10	FAL23734
MDD	POL-FR-03-23-230-SU	3 Prong	2.3	230	10	FAL23530
MDD	POL-FR-04-23-230-SU	4 Prong	2.3	230	10	FAL23778
MDD	POL-FR-05-23-230-SU	5 Prong	2.3	230	10	FAL23643
Pentax	FG-23-230-A1	Alligator	2.8	230	10	
Pentax	FG-23-230-A2	Rat Tooth	2.8	230	10	

# ENDO DECONTAMINATION



# Single-Ended Cleaning Brushes

## Cleaning & Disinfection



- Robust durable cleaning brush provides a smooth channel pass and confidence in removal of soil from inside the working channel
- Uniformed bristle density offers a consistent lumen contact to remove soil from channel


Supplier	Product Code	Operating Channel Diameter (mm)	Catheter		Brush (mm)		QOM	FAL Code
			Length & Colour (cm)	Diameter (mm)	Single / Double / Triple	Length		
<b>Adult: gastroscope, colonoscope, duodenoscope, enteroscope and cystoscope and EUS endoscopes</b>								
Vytil	CJ-ESB-020-06	2.8 to 5.0	20 Light Blue	1.8	—●	20	6	50 FAL24651
Vytil	CJ-ESB-230-06	2.8 to 5.0	230 Blue	1.8	—●	20	6	50 FAL23706
MDD	SCBS-23-20-06-230	1.7	230 Yellow	1.7	—●	20	6	100 FAL23764

# Double-Ended Cleaning Brushes



Supplier	Product Code	Operating Channel Diameter (mm)	Catheter		Brush (mm)		QOM	FAL Code
			Length & Colour (cm)	Diameter (mm)	Single / Double / Triple	Length		
Vytil	CJ-EDB-100-02	1.2 to 1.5	100 White	1	●—●	10 / 10	2	50 FAL24656
Vytil	CJ-EDB-120-02	1.2 to 1.5	120 White	1	●—●	10 / 10	2	50 FAL24649
Vytil	CJ-EDB-180-02	1.2 to 1.5	180 White	1	●—●	10 / 10	2	50 FAL24635
Vytil	CJ-EDB-230-02	1.2 to 1.5	230 White	1	●—●	10 / 10	2	50 FAL24631
Vytil	CJ-EDB-120-03	2.0 to 2.5	120 Yellow	1.6	●—●	20 / 20	3	50 FAL24639
Vytil	CJ-EDB-230-03	2.0 to 2.5	230 Yellow	1.6	●—●	20 / 20	3	50 FAL24647
Vytil	CJ-EDB-300-03	2.0 to 2.5	300 Yellow	1.6	●—●	20 / 20	3	50
Vytil	CJ-EDB-120-06	2.8 to 5.0	120 Blue	1.8	●—●	20	6	50 FAL24638
Vytil	CJ-EDB-230-06	2.8 to 5.0	230 Blue	1.8	●—●	20 / 20	6	50 FAL24640
Vytil	CJ-EDB-300-06	2.8 to 5.0	300 Blue	1.8	●—●	20 / 20	6	50
MDD	DCBS-08-10-15-120	1.0 to 1.5	120 Blue	0.8	●—●	10 / 10	1.5	100 FAL23792
MDD	DCBS-08-10-15-180	1.0 to 1.5	180 Blue	0.8	●—●	10 / 10	1.5	100 FAL23899
MDD	DCBS-08-10-15-230	1.0 to 1.5	230 Blue	0.8	●—●	10 / 10	1.5	100 FAL23745
MDD	DCBS-15-15-02-120	1.2 to 1.8	120 Yellow	1	●—●	15 / 15	2	100 FAL23906
MDD	DCBS-15-15-06-120	1.2 to 1.8	120 Yellow	1	●—●	15 / 15	6	100 FAL23539
MDD	DCBS-15-15-02-230	1.2 to 1.8	230 Blue	1	●—●	15 / 15	2	100 FAL23779
MDD	DCBS-18-20-03-120	2.0 to 2.8	120 Green	1.5	●—●	20 / 20	3	100 FAL23877
MDD	DCBS-18-20-03-180	2.0 to 2.8	180 Green	1.5	●—●	20 / 20	3	100 FAL23499
MDD	DCBS-18-20-03-230	2.0 to 2.8	230 Green	1.5	●—●	20 / 20	3	100 FAL23579
MDD	DCBS-23-20-06-230	2.8 to 4.2	230 Blue	1.7	●—●	20 / 20	6	100 FAL23648
MDD	DCBS-23-20-06-300	2.8 to 4.2	300 Blue	1.7	●—●	20 / 20	6	100 FAL23614
MDD	DCBS-23-20-05-180	2.8 to 4.2	180 Yellow	1.8	●—●	20 / 20	5	100 FAL23765
MDD	DCBS-23-20-05-230	2.8 to 4.2	230 Yellow	1.8	●—●	20 / 20	5	100 FAL23543
MDD	DCBS-23-20-05-180-X22		180 Green	2.2	●—●	20 / 20	6	100 FAL23555
MDD	DCBS-23-20-06-230-X22		240 Green	1.8	●—●	20 / 20	6	100 FAL23586
MDD	DCBS-23-20-06-240		240 Blue	1.8	●—●	20 / 20	6	100 FAL23597

# Multi-Point Cleaning Brushes



- Ball heads on the brush ends and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Ideal brush head size and density
- Brush-head diameter allows for maximum mechanical brushing of the endoscope channel
- Bend resistant thanks for first-class catheter rigidity

Supplier	Product Code	Operating Channel Diameter (mm)	Catheter		Brush (mm)			QOM	FAL Code
			Length & Colour (cm)	Diameter (mm)	Single / Double / Triple	Length	Diameter		
MDD	MCBS-23-20-06-230	2.8 to 4.2	230 Blue	1.7	●—●—●	20 / 20 / 20	6	100	FAL23861
MDD	MCBS-23-20-06-260	2.8 to 4.2	260 Blue	1.7	●—●—●	20 / 20 / 20	5	100	FAL23771

# Triple-Point Cleaning Brushes



- The triple brush head cleaning device covers 3x more contact surface area in comparison to a standard cleaning brush
- Has both distal & proximal valve brush attachments
- Atraumatic tip to prevent channel scarring
- Reinforced catheter provides both flexibility while maintaining somewhat stiff to ensure ease of passage

Supplier	Product Code	Operating Channel Diameter (mm)	Catheter		Brush (mm)			QOM	FAL Code
			Length & Colour (cm)	Diameter (mm)	Single / Double / Triple	Length	Diameter		
MDD	TCBS-23-20-06-230	2.8 to 4.2	230 Blue	1.7	●—●—●	20 / 20 / 125	6	100	
Vytil	NLS/ETB-230-06505	2.6 to 5.0	230 Blue	1.8	●—●—●	20 / 20 / 125	6	100	FAL24637

# Sweeper Brushes



- Sweeper brush gives maximum contact with the internal channels of the endoscope, ensuring debris and bio-burden is swept out in a single pass
- Ball heads on the brush ends to prevent brush shed, and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon with brush-head diameter allowing for maximum mechanical brushing of the endoscope channel
- Full 360-degree circumferential lumen contact cleaning increases cleaning efficacy
- Easy one direction pass through ensures contents are pulled through the channel and not reintroduced back into the lumen.

Supplier	Product Code	Operating Channel Diameter (mm)	Catheter		Brush (mm)		Sweeper (mm)		QOM	FAL Code
			Length & Colour (cm)	Diameter (mm)	Length	Diameter	Length	Diameter		
MDD	SCBS-18-20-3-300-3-120	1.8 to 2.8	120 Green	1.5	20	3	300	3	100	FAL23500
MDD	SCBS-18-20-3-300-3-160	1.8 to 2.8	160 Orange	1.5	20	3	300	3	100	FAL23580
MDD	SCBS-23-20-5-300-5-160	2.8 to 4.2	160 Blue	1.7	20	5	300	5	100	FAL23719
MDD	SCBS-23-20-5-300-5-230	2.8 to 4.2	230 Yellow	1.7	20	5	300	5	100	FAL23702

# Combination Cleaning Brushes



DVBS-23-40-10-20-06-230  
Without Handle




DVEC-23-40-10-15-05-230  
With Handle

- Ball heads on the brush ends to prevent brush shed, and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon with brush-head diameter allowing for maximum mechanical brushing of the endoscope channel

Supplier	Product Code	Operating Channel Diameter (mm)	Catheter		Valve Brush (mm)		Channel Brush (mm)		QOM	FAL Code
			Length & Colour (cm)	Diameter (mm)	Length	Diameter	Length	Diameter		
MDD	DVBS-23-40-10-20-06-230	2.8 to 4.2	230 Light Blue	1.7	40	10	20	6	100	FAL23854
MDD	DVEC-23-40-10-15-05-230	2.8 to 4.2	230 White	1.7	40	10	15	5	100	FAL23645
MDD	DVBS-23-40-11-20-05-230	2.8 to 4.2	230 Blue	1.7	40	11	20	5	100	FAL23693
MDD	DVBS-23-40-11-20-06-230	2.8 to 4.2	230 Blue	1.7	40	11	20	6	100	FAL23884

# X-MAS Cleaning Brushes



DVBS-23-20-07/4.5-20-05-230

- Christmas tree shaped valve brushes to enable cleaning in hard to reach places of the scope
- Sweeper brush gives maximum contact with the internal channels of the endoscope, ensuring debris and bioburden is swept out in a single pass
- Ball heads on the brush ends to prevent brush shed, and rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon with brush-head diameter allowing for maximum mechanical brushing of the endoscope channel

## Two-Sided X-MAS Tree Brush with Sweeper

Supplier	Product Code	Catheter		Valve Brush (mm)		Sweeper Brush (mm)		Channel Brush (mm)		QOM	FAL Code
		Length & Colour (cm)	Diameter (mm)	Length	Diameter	Length	Diameter	Length	Diameter		
<b>2.8 to 4.2 Operating Channel Diameter (mm)</b>											
MDD	DVBS-23-20-07/4.5-20-05-230	Blue 230	1.7	20	4.5 - 7	100	5	10	3 - 6	100	FAL23497
MDD	DVBS-23-20-07/4.5-20-05-260	Blue 260	1.7	20	4.5 - 7	100	5	10	3 - 6	100	FAL23549

## Two-Sided X-MAS Tree Brush

Supplier	Product Code	Operating Channel Diameter (mm)	Catheter		Valve Brush (mm)		Channel Brush (mm)		QOM	FAL Code
			Length & Colour (cm)	Diameter (mm)	Length	Diameter	Length	Diameter		
MDD	DCBS-23-20-07/4.5-20-05-230	2.8 to 4.2	Blue 230	1.7	20	4.5 - 7	10	3 - 6	100	FAL23830

# Valve Brushes



DVS-30-35-10-20-05-15



SVS-30-35-11



CJ-GDB-11-05



CJ-GDB-11-05

- Rounded brush tips for gentle cleaning
- Brushes made of soft, durable nylon
- Rounded handle

Supplier	Product Code	Catheter	Valve Brush (mm)		Channel Brush (mm)		QOM	FAL Code
		Length & Colour (cm)	Length	Diameter	Length	Diameter		
MDD	DVS-30-35-10-20-05-15	15 Blue	35	11	20	5	200	FAL23859
MDD	SVS-30-35-11	9.5 Green	35	11	-	-	200	FAL23564
MDD	SES-30-2-15	8 Red	-	-	15	2	200	FAL23851
Vytil	CJ-GDB-11-05	20 White	32 / 20	11 / 5	-	-	50	FAL24644
Vytil	CJ-GSB-02	9 White	10	2	-	-	50	FAL23886



- U-shape brush for endoscope sheath
- Compatible with 22 to 28 mm endoscope sheath

Supplier	Product Code	Axis		Brush (mm)		QOM	FAL Code
		Length (cm)	Diameter (mm)	Length	Diameter		
Vytil	NLS/Ubrush	150	4	80	25	25	

Supplier	Product Code	Operating channel diameter (mm)	Catheter		Brush (mm)		QOM	FAL Code	
			Length & Colour (cm)	Diameter (mm)	Simple, Double or Triple	Length			Diameter
<b>1 Cleaning brush and 1 Valve cleaning brush</b>									
Vytil	CJ-KEDB-100-02	1.2 to 1.5	100 - White	1.0	●—●	10 / 10	2	40	FAL24659
Vytil	CJ-KEDB-120-02	1.2 to 1.5	120 - White	1.0	●—●	10 / 10	2	40	FAL24630
Vytil	CJ-KEDB-180-02	1.2 to 1.5	180 - White	1.0	●—●	10 / 10	2	40	
Vytil	CJ-KEDB-120-03	2.0 to 2.5	120 - Yellow	1.6	●—●	20 / 20	3	40	FAL24641
Vytil	CJ-KEDB-230-03	2.0 to 2.5	230 - Yellow	1.6	●—●	20 / 20	3	40	FAL24655
Vytil	CJ-KEDB-300-03	2.0 to 2.5	300 - Yellow	1.6	●—●	20 / 20	3	40	
MDD	DCB-SET-23-20-07/ 4.5-20-05-230	2.0 to 2.5	230 - Blue	1.7	●—●	20 / 20	3 / 6	100	FAL23662
MDD	DCBS-SET-23-20-3/ 6-10-4/7-230		230 - Green	1.65	●—●	20 / 20	3 / 6	100	FAL23802
MDD	DCBS- SET-08-10-15-165		165 - Blue	0.8	●—●	10 / 10	1.5	100	FAL23621
MDD	DCBS- SET-23-20-05-240		240 - Pink	1.65	●—●	10 / 10	1.5	100	FAL23625
MDD	DCBS-SET-180 -20-05-35-10		180 - Yellow	1.8	●—●	20 / 20	5	100	FAL23903
MDD	DCBS- SET-20-05-35-10		230 - Yellow	1.8	●—●	20 / 20	5	100	FAL23795

**Adult: Sets for gastroscope, colonoscope, duodenoscope, enteroscope, cystoscope and EUS scope**

Vytil	CJ-KEDB-120-06	2.8 to 5.0	120 - Blue	1.8	●—●	20 / 20	6	40	FAL24654
Vytil	CJ-KESB-230-06	2.8 to 5.0	230 - Blue	1.8	—●	20	6	40	FAL24645
Vytil	CJ-KEDB-230-06	2.8 to 5.0	230 - Blue	1.8	●—●	20 / 20	6	40	FAL24658
Vytil	CJ-KEDB-300-06	2.8 to 5.0	300 - Blue	1.8	●—●	20 / 20	6	40	

**Triple cleaning brush including a sweeper brush**

Vytil	NLS/KETB-230-06S05	2.8 to 5.0	230 - Blue	1.8	●—●—●	20 / 20 / 125	6	40	FAL24632
MDD	MCBS- SET-23-20-06-230	2.8 to 4.2	230 - Blue	1.7	●—●—●	20 / 20 / 20	6	100	FAL23883
MDD	MCBS- SET-23-20-06-260	2.8 to 4.2	260 - Blue	1.7	●—●—●	20 / 20 / 20	6	100	FAL23654
MDD	SCBS- SET-20-05-300-5-230	2.8 to 5.0	230 - Blue		●—●—●	20 / 20 / 300	5	100	FAL24558

**2 Cleaning brushes and 1 Valve cleaning brush**

Vytil	CJ-KE-230-06- 230-02	1.2 to 1.5	230 - White	1.0	●—●	10 / 10	2	40	FAL24650
		2.8 to 5.0	230 - Blue	1.8	●—●	20 / 20	6		
		/	20 - White	/	●=●	32 / 20	11 / 5		
Vytil	CJ-KE-120-03- 020-06	2.0 to 2.5	120 - Yellow	1.6	●—●	20 / 20	3	40	FAL24642
		2.8 to 5.0	20 - Light Blue	1.8	—●	20	6		
		/	20 - White	/	●=●	32 / 20	11 / 5		





**Protective Bag:**

- To protect the endoscope or the transport tray
- Saves space when rolled up
- Self-closing, flexible elastic band

**CC-Transport Bag for Turning Around:**

- “Clean” and “Contaminated” marking on transparent polyethylene film (PE)
- Saves space when rolled up
- Self-closing, flexible elastic band

Supplier	Product Code	Dimensions (cm)		QOM	FAL Code
		Length	Width		
MDD	STB-01-SET	70	58	25	FAL23473
MDD	STB-02-SET	95	80	25	FAL23694



Our range of distal tip protectors help to protect endoscopes during transportation & storage from both potentially costly repairs & downtime

- Simple & easy to use design
- Reduce chance of damage and costly repairs to endoscope
- Disposable & compatible for all sizes & brands of endoscopes

Supplier	Product Code	Length (mm)	Suitable For Instruments with a Diameter of (mm)	Material	Sterile	QOM	FAL Code
MDD	DTP-01	50	7 - 15	Mesh	Yes	100	FAL23504
MDD	DTP-02	50	2 - 10	Mesh	Yes	100	FAL23761
MDD	EP-05-01	140	4 - 24	Sponge	No	50	FAL23685
MDD	EP-05-01-EO	140	4 - 24	Sponge	Yes	50	FAL23696
MDD	EP-14-01	140	9 - 24	Sponge	No	50	FAL23731
MDD	EP-14-01-EO	140	9 - 24	Sponge	Yes	50	FAL23560



- Optimal protection of endoscopes during transport
- Secure support of endoscopes during storage
- Reliable cleaning of endoscopes

Supplier	Product Code	Length (mm)	Outer Diameter (mm)	Sterile	QOM	FAL Code
MDD	EP-05-14-CL	140	37	No	50	FAL23615
MDD	EP-05-14-CL-EO	140	37	Yes	50	FAL23887
MDD	EP-14-01-EO	140	37	Yes	50	FAL23560

# GI STENTING






AN/ESS-20080-F



AN/ESA-18090-F


- Comprehensive range of Oesophageal self-expanding metal stents provides clinicians with stent design options to suit specific indications
- High density Nitinol Wire ensures strong residual force whilst in situ in the lumen
- Anti-migration features help mitigate migration risks of the stent in situ
- 11 Tantalum Radiopaque markers provide the end user with excellent visibility

Supplier	Product Code	Description	Total Length (mm)	Diameter Middle / Ends (mm)	Shape	QOM	FAL Code
Vytil	AN/ESA-18090-F	Covered	9	18 / 24	Single Spherical	1	FAL23805
Vytil	AN/ESA-18110-F	Covered	11	18 / 24	Single Spherical	1	FAL23545
Vytil	AN/ESS-20080-F	Covered	8	20 / 26	Double Cup	1	FAL23806
Vytil	AN/ESS-20110-F	Covered	11	20 / 26	Double Cup	1	FAL23807
Vytil	AN/ESS-20140-F	Covered	14	20 / 26	Double Cup	1	FAL23684
Vytil	AN/ESS-20100-P	Partially Covered	10	20 / 26	Double Cup	1	FAL23814
Vytil	AN/ESS-20150-P	Partially Covered	15	20 / 26	Double Cup	1	FAL23740




- Application in palliative treatment of colorectal strictures caused by malignant tumours.
- Dual layered mesh design with very small cavities to prevent tissue ingrowth
- Very low returning force provides excellent conformity to patient anatomy whilst significantly reducing stent migration and risk of micro perforation
- Atraumatic ends reduce ulcerations and discomfort
- Distal lasso to aid accurate repositioning and retrieval
- Reloadable delivery device; up to 70% of the stent can be deployed before re sheathing
- TTS 230cm delivery system

Supplier	Product Code	Description	Total Length (mm)	Introducer Diameter (Fr)	Shape	QOM	FAL Code
Vytil	AN/IS-22060-N	Uncovered	60	10	Double Cup	1	FAL24551
Vytil	AN/IS-22090-N	Uncovered	90	10	Double Cup	1	FAL24556
Vytil	AN/IS-22120-N	Uncovered	120	10	Double Cup	1	FAL24567



- Round tip for alpha loop forming
- Black and white stripes on 30 cm to display wire movements
- Radio-opaque hydrophilic distal end on 50 - 55 mm
- Nitinol core for good stiffness / push compromise

Supplier	Product Code	Description	Flexibility	Length (cm)	Diameter (mm / inch)	QOM	FAL Code
Vytil	AN/CGW-S3526	Straight guide wire	Standard	260	0.89 / 0.0035"	10	FRD6934
Vytil	AN/CGW-S3526S		Stiff	260	0.89 / 0.0035"	10	FRD6935
Vytil	AN/CGW-S3545		Standard	450	0.89 / 0.0035"	10	FRD6938
Vytil	AN/CGW-S2526S		Stiff	260	0.89 / 0.0025"	10	
Vytil	AN/CGW-S3545S		Stiff	450	0.89 / 0.0035"	10	FRD6936
Vytil	AN/CGW-S2545S		Stiff	450	0.89 / 0.0035"	10	FRD6941
Vytil	AN/CGW-A2526S	Angled guide wire	Stiff	260	0.89 / 0.0025"	10	
Vytil	AN/CGW-A3526		Stiff	260	0.89 / 0.0035"	10	FRD6939
Vytil	AN/CGW-A3545		Standard	450	0.89 / 0.0035"	10	FRD6937
Vytil	AN/CGW-A2545S		Stiff	450	0.63 / 0.0025"	10	FRD6940



### Savary with nitinol core

- Super rigid nitinol guidewire, flexible distal end of 10 cm
- Good push/flexibility compromise
- To use for a stent placement or an esophagus dilation

### Savary type

- Super rigid stainless steel guidewire, flexible distal end of 10 cm
- To use for a stent placement or an esophagus dilation

### Super stiff

- Super stiff and coated metallic guidewire, flexible distal end of 9 cm
- To use for a stent placement or an esophagus, colonic, enteral dilation

Supplier	Product Code	Description	Flexibility	Length (cm)	Diameter (mm / inch)	QOM	FAL Code
Vytil	47012263	Savary with nitinol core guide wire		230	0.035"	10	
Vytil	43002341	Savary type guide wire		230	0.035"	10	
Vytil	46004758	Super stiff guidewire	Stiff	460	0.035"	10	

# CAPSULE ENDOSCOPY




# Small Bowel Capsule Endoscopy



AKES-11SW



SensorBelt



ESNavi

The NaviCam System with ProScan, based on a deep-learning algorithm has a lesion detection sensitivity of 99.9% and a 5.9 minute average read time per patient.

The NaviCam range combines exceptional technology and optics with innovative integration of Artificial Intelligence for assisted diagnosis.

NaviCam comprises the first Capsule Endoscopy system with the addition of ProScan, a deep-learning CNN platform able to distinguish and detect pathologies across over 12 variations. With data published, and presented on the cover page, of Gastroenterology 2019, the AI model has a lesion detection sensitivity of 99.9% and a reduced read time to under 6 minutes.

Supplier	Product Code	Description	QOM	FAL Code
Anx Robotica	AKS-1	Capsule Locator	1	FAL18995
Anx Robotica	PackingCase	Metal Packing Case	1	FAL18992
Anx Robotica	SensorBelt	Sensor Belt Gen 1	1	FAL19001
Anx Robotica	AKES-11SK	NaviCam™ SB Capsule Gen 2	10	
Anx Robotica	AKES-11SW	NaviCam™ SB Capsule Gen 1	10	FAL18996
Anx Robotica	ESNavi	ESNavi Software	1	FAL19002
Anx Robotica	ProScan	NaviCam™ ProScan Full Licence	1	FAL18994
Anx Robotica	AKR-2	Data Recorder Belt	1	
Anx Robotica	NS-1	Digital Recorder Host	1	FAL19000
Anx Robotica	NS-1/AKRK-2	NaviCam™ SB Capsule Endoscopy System With Data Recorder, Locator & Software	1	

# Stomach Capsule Endoscopy



- The first 360° robotic-controlled stomach capsule-revolutionising capsule endoscopy
- The risk of perforation and cross-contamination is minimized with MCCE
- NaviCam can navigate a 180 degree turn in order to view the oesophageal sphincter
- MCCE can view the stomach in its natural state rather than the sedated invasive state, which is important for the diagnosis of motility conditions such as gastroparesis
- Because NaviCam Stomach Procedures do not require sedation, patients can resume normal activities following the procedure. NaviCam also provides an alternative to those patients in which sedation is contraindicated due to co-morbidities, or health risk.

Supplier	Product Code	Description	QOM	FAL Code
Anx Robotica	AKS-1	Capsule Locator	1	FAL18995
Anx Robotica	PackingCase	Metal Packing Case	1	FAL18992
Anx Robotica	AKEM-11SW	NaviCam™ Stomach Capsule	10	FAL18993
Anx Robotica	ESNavi	ESNavi Software	1	FAL19002
Anx Robotica	ProScan	NaviCam™ ProScan Full Licence	1	FAL18994
Anx Robotica	AKR-2	Data Recorder Belt	1	
Anx Robotica	NC-1	NaviCam™ Stomach Capsule System With Controller, Data Recorder, Locator & Software	1	FAL18998

MEDI TECH







The minimally invasive endoluminal approach has evolved to become a mainstream treatment option for obesity. Particularly, endoscopic remodelling of the stomach and restricting the gastric lumen by placing multiple, internal, full-thickness suture plication has been shown to affect hunger and appetite and to induce effective and sustained weight-loss. The modified POSE procedure (POSE 2.0) aims to narrow and shorten the stomach significantly in comparison with its original size. Involving the use of the incisionless operating platform, POSE 2.0 has 4 distinct components:

- Transport
- g-Lix
- g-Prox EZ
- g-Cath EZ

Supplier	Product Code	Description	QOM	FAL Code
USGI	118180	Transport Retroflex	1	
USGI	208100	g-Cath EZ Delivery Catheter	6	
USGI	208200	g-Prox EZ 16	1	
USGI	208220	g-Prox EZ 33	1	
USGI	301010	g-Lix Flexible Tissue Anchor Catheter	3	
USGI	500100	Transport-R (Right Angle)	1	
USGI	POSE2.0 Procedure Kit	POSE2.0 Procedure Kit Contains: 500100, 208220, 301010, 208100	1	



- Demonstrated to be safe and effective
- Precision delivery of non-thermal ablation which targets tumour tissue whilst preserving surrounding healthy tissue structures
- ePORE and EndoVE are repeatable and leave standard treatment options available
- Delivered under light sedation, day case, cost effective
- Simple standard connection for endoscopes between 9mm and 10mm
- Connects to a vacuum which draws the tumour tissue into the chamber this can be viewed via the window at the head of the electrode
- Ease of use

Supplier	Product Code	Description	QOM	FAL Code
Mirai Medical	MS-18-07	EndoVE High Frequency Extension Lead	1	
Mirai Medical	MS-19-05	EndoVE Low Frequency Extension Lead	1	
Mirai Medical	Epore MMHT1	ePORE Generator MMHT1	1	
Mirai Medical	EndoVE 700209	EndoVE 700209 Electroporation Probe	1	



Taking feedback from experienced clinicians across Europe our CUTIS skin electrode has been designed and engineered with both the clinician and patient in mind.

The retractable head offers unique benefits:

- Needles are fully supported as they enter and penetrate tissue, avoiding needle bending, convergence
- Patients treated under LA do not see the needles, relieving potential anxiety
- Safe - no need for recapping

CUTIS AND ePORE® ADVANTAGES

- High frequency electroporation with no muscular contraction
- Ability to be delivered under LA in an outpatient setting
- Precision delivery of non-thermal ablation
- Preservation of healthy tissue
- CUTIS and ePORE® are repeatable and leave standard treatment options available
- Significant cost savings

Supplier	Product Code	Description	QOM	FAL Code
Mirai Medical	MS-19-03	CUTIS High Frequency Extension Lead	1	
Mirai Medical	MS-19-04	CUTIS Low Frequency Extension Lead	1	
Mirai Medical	Epore MMHT1	ePORe Generator MMHT1	1	
Mirai Medical	Cutis 700245	CUTIS 700245 Probe	5	

CADDIE is Odin Vision's cutting-edge Artificial Intelligence system for a new era of endoscopy. It supports endoscopists to detect and characterise colorectal polyps in real-time during colonoscopy procedures:

AI to improve colonoscopy:

- AI Polyp Detection
- AI Polyp Characterisation
- AI Bowel Cleanliness (Visible Mucosa)
- AI Caecum (AO Confirmation)
- Cloud Technology
- Vendor Neutrality
- Easy to Use

Supplier	Product Code	Description	QOM	FAL Code
Odin Vision	OdinACT001	Set Up and Activation Cost for Consignment	1	
Odin Vision	OdinACT002	Bronze Subscription Set Up and Activation Cost	1	
Odin Vision	OdinACT003	Silver Subscription Set Up and Activation Cost	1	
Odin Vision	OdinACT004	Gold Subscription Set Up and Activation Cost	1	
Odin Vision	OdinPAG001	Bronze Pay as You Go Detection Only Micro up to 250 procedures, pay per procedure	1	
Odin Vision	OdinPAG002	Silver Pay as You Go Detection Only Small up to 500 procedures, pay per procedure	1	
Odin Vision	OdinPAG003	Gold Pay as You Go Detection Only Medium up to 1500 procedures, pay per procedure	1	
Odin Vision	OdinPAG004	Platinum Pay as You Go Detection Only Large up to 3000 procedures, pay per procedure	1	
Odin Vision	OdinPAG005	Bronze Pay as You Go Detection and Diagnosis Micro up to 250 procedures, pay per procedure	1	
Odin Vision	OdinPAG006	Silver Pay as You Go Detection and Diagnosis Small up to 500 procedures, pay per procedure	1	
Odin Vision	OdinPAG007	Gold Pay as You Go Detection and Diagnosis Medium up to 1500 procedures, pay per procedure	1	
Odin Vision	OdinPAG008	Platinum Pay as You Go Detection and Diagnosis Large up to 3000 procedures, pay per procedure	1	
Odin Vision	OdinPC01	Hardware for Consignment	1	
Odin Vision	OdinPC02	Bronze Subscription Hardware	1	
Odin Vision	OdinPC03	Silver Subscription Hardware	1	
Odin Vision	OdinPC04	Gold Subscription Hardware	1	
Odin Vision	OdinSUB001	Bronze Subscription Detection Only up to 500 procedures 2 rooms installed, pay per procedure	1	
Odin Vision	OdinSUB002	Silver Subscription Detection Only up to 1500 procedures 3 rooms installed, pay per procedure	1	
Odin Vision	OdinSUB003	Gold Subscription Detection Only up to 3000 procedures all rooms installed, pay per procedure	1	
Odin Vision	OdinSUB004	Bronze Subscription Detection and Diagnosis up to 500 procedures 2 rooms installed, pay per procedure	1	
Odin Vision	OdinSUB005	Silver Subscription Detection and Diagnosis up to 1500 procedures 3 rooms installed, pay per procedure	1	
Odin Vision	OdinSUB006	Gold Subscription Detection and Diagnosis up to 3000 procedures all rooms installed, pay per procedure	1	

# CLINICAL TRAINING





The T(ractus) D(igestives) Box (TD-Box) makes hard- or showcase training easy and professional. For instance, it brings medical training with (animal) models to the next level. With a compact training box equipped with dark and transparent covers, convenient storage facilities and is compatible with all major electrocautery brands.

Our compact training box makes hard- or showcase training easy and flexible. Portable, compatible with all major electrocautery brands and equipped with storage facilities, this made-for-purpose kit perfectly fits the complicated requirements of modern training curriculums.

Supplier	Product Code	Description	QOM	FAL Code
Medical Training Tools	100.01.01	TD Box Mobile	1	
Medical Training Tools	200.01.01	Upper GI Model	1	
Medical Training Tools	200.01.02	Transparent Upper GI Model	1	
Medical Training Tools	300.01.01	Cleaning Kit	1	
Medical Training Tools	400.01.01	Foreign Body Removal Kit	1	
Medical Training Tools	700.01.01	Biliary Tree Kit	1	
Medical Training Tools	800.01.01	Self-Sticking Velcro	1	
Medical Training Tools	700.01.02	Biliary Tree Kit (Non Transparant)	1	

# WATCHMAN MANAGED INVENTORY



# Watchman Managed Inventory

Duomed UK are driven to making our customer's life easier through dedicated service and support. Our Watchman Managed Inventory can help you reduce the cost of holding unnecessary inventory and streamline stock management within your hospital. In addition to working hand-in-hand with your procurement and clinical teams, we seamlessly ensure you have what you need, when you need it.

Following inventory levels being forecast and agreed upon, stock is supplied on a consignment basis. Ownership of stock remains with Duomed UK, with responsibility of safe storage moving to the end user. Devices are only invoiced once used by the end user and invoiced accordingly. The key benefits to our Watchman service include:

- Freeing up liquid asset being held on shelves, allowing investment in other areas of your department/hospital
- Mitigation of short date risk, as Duomed UK remain responsible for stock management, checks and rotation
- Reduced responsibility and time of nurse/staff support to check, and order required stock levels, permitting for time to focus on other key task areas
- Reduction in overall purchase orders generated, with a subsequent cost efficiency. Average cost to generate one purchase order is circa £75

## Benefits of Partnership with Duomed UK

Time is something we are always short of and is difficult to place a financial value on. Watchman Managed Inventory offers the ability to save hospital working hours lost due to stock management and lengthy invoicing runs.

Our commercial team will work with you utilising a monthly tracker, allowing clear visibility to spend management and stock control. This eliminates the need for the hospital staff to carry out large numbers of stock checks and reduces the number of purchase orders that are needed. Consolidated invoicing through one supplier also reduces the number of payment processes required to efficiently service an Endoscopy department.

### Mitigate risk

Rationalising the number of suppliers enables familiarity with device brands and functionality. Avoiding variance promotes confidence and competence in supplier and device function, reducing risks associated with unfamiliarity.

### Better patient outcomes

Confidence and competence in brand and devices supplied instils a level of experience and expertise. When used consistently through reproducible application, this builds best practice and improves the likelihood of better patient outcomes.

### Access to innovation

A partnership with Duomed UK enhances the ability to access innovation and new product launches. To access new products, you simply need to contact your Endoscopy Portfolio Manager and ask for the product to be trialled.

### Improved referral to treatment times

Through a combination of investment in new niche innovation and best practice through familiarity, RTT times can be improved through earlier diagnosis and better patient pathway access. Efficiency drives better outcomes and spend efficiencies through all department areas.