

Hot Biopsy Forceps



- This forcep is used for small polyps removal
- Oval and fenestred jaws made of surgical stainless steel
- PTFE coated catheter
- Coagulation is achieved with opened or closed jaws

Supplier	Product Code	Forceps		Catheter		QOM	FAL Code
		Jaws	Spike	Length (cm)	Diameter (mm)		
Vytil	CJ-PC-23-230	= standard	No	230	2.3	10	FAL23688
Vytil	CJ-PCA-23-230	= standard	Yes	230	2.3	10	FAL23583
Vytil	NLS/HotBFO-24-230	= standard	No	230	2.4	10	