


Tracheo-Bronchial Biopsy Forceps



Alligator jaw with spike



Standard jaw without spike

- Graduated catheter to easily identify the biopsy forceps progression within the operating channel
- All forceps are delivered with a sample retriever to reduce blood exposure risk

Supplier	Product Code	Forceps			Catheter		QOM	FAL Code
		Jaws	Spike	PTFE coated	Length (cm)	Diameter (mm)		
Vytil	AN/AMHBFA	= standard	No	Yes	110	1.8	10	
Vytil	AN/AMHBFA-WC	= standard	No	No	110	1.8	10	
Vytil	AN/AMHBFC	≡ alligator	No	Yes	110	1.8	10	
Vytil	AN/AMHBFC-WC	≡ alligator	No	No	110	1.8	10	
Vytil	CJ-PPT-18-110-A	= standard	Yes	Yes	110	1.8	10	FAL24646