



Spiral snare

- Spiral snare has additional spiral wire around the snare to prevent slipping on the mucosa
- Spiral snare has an ultra stiff wire
- Also a standard hot rotatable snare is available

Supplier	Product Code	Snare Width (mm)	Catheter		QOM	FAL Code
			Diameter (mm)	Length (cm)		
Vytil	CJ-ADR-23-230-HCO06	6	2.3	230	10	
Vytil	VTL/SN-24-240-HCO10	10	2.3	230		FAL23678
Vytil	CJ-ADR2-23-230-NG20	20 x 25	2.3	230		FAL23619
Vytil	CJ-ADR2-23-230-NG25	25 x 35	2.3	230		FAL23782
Vytil	CJ-ADR2-23-230-NG30	30 x 45	2.3	230		FAL23889