



## Superior Navigation

- Inflatable two-balloon design stabilises the endoscope and the nearby anatomy
- Intraoperative conduit enables safe and efficient access, and repeated removal of tools and tissue, without losing position

## Enhanced Visualisation

- Aft balloon supplies gentle endoluminal positioning, helping to reveal tissue planes and expose hidden polyps within folds and flexures.
- EZ Glide hydrophilic coating provides a clearer, more visible therapeutic zone without messy lubricants

## Stability & Control

- Colon shortening allows for safe and gentle straightening of tortuous anatomy, achieving faster access and fewer aborted procedures
- Extendable fore balloon creates a stable therapeutic zone, flattening the surgical site for submucosal injection and precise tissue dissection
- Integrated suture system supplies controlled, dynamic tissue manipulation and tension during procedures

Supplier	Product Code	Length (cm)	QOM	FAL Code
Lumendi	D-1001-5	130	5	
Lumendi	D-1002-5	103	5	
Lumendi	D-1001	130	1	
Lumendi	D-1002	103	1	