

Dilation Balloons



Dilation



- Available for oesophagus, cardia, pylorus, colon and biliary duct.
- All balloons are TTS compatible except achalasia balloons
- One size fits all
- 2 radiopaque marks
- Compatible with guidewires of 0.035 "maximum

Supplier	Product Code	Balloon		Catheter		Pressure (ATM) Recommended / Burst	Comments	QOM	FAL Code
		Diameter (mm)	Length (mm)	Diameter (mm/fr)	Length (cm)				
Dilation balloon single diameter									
Vytil	CJ-JHY-BD-06-40-250	6	40	2.4 / 7	250	8 - 10			FAL24573
Vytil	CJ-JHY-BD-08-40-250	8	40	2.4 / 7	250	8 - 10			FAL24593
Vytil	CJ-JHY-BD-10-40-250	10	40	2.4 / 7	250	8 - 10	Operating channel diameter: ≥ 2.8 mm	1	FAL24577
Vytil	CJ-JHY-BD-12-60-250	12	60	2.4 / 7	250	5 - 6			FAL24568
Vytil	CJ-JHY-BD-14-60-250	14	60	2.4 / 7	250	5 - 6	Guide wire: 0.035"		FAL24557
Vytil	CJ-JHY-BD-16-60-250	16	60	2.4 / 7	250	5 - 6		FAL24553	
Vytil	CJ-JHY-BD-16-80-250	16	80	2.4 / 7	250	5 - 6			FAL24550
Vytil	CJ-JHY-BD-18-80-250	18	80	2.4 / 7	250	4 - 5			FAL24580
Vytil	CJ-JHY-BD-20-80-250	20	80	2.4 / 7	250	4 - 5			FAL24552
Achalasia balloon single diameter									
Vytil	CJ-JHY-BD-30-80-90	30	80	6.7 / 20	90	2 - 3	OTW procedure guide wire: 0.035"	1	FAL24596
Vytil	CJ-JHY-BD-35-80-90	35	80	6.7 / 20	90	2 - 3			FAL24562
Vytil	CJ-JHY-BD-40-80-90	40	80	6.7 / 20	90	2 - 3			FAL24565