


# Monofilament Snares



CJ-ADR2-23-230-O25B



CJ-ADR-23-230-MFO-O20

- Stiff single element monofilament wire provides maximum snare rigidity for challenging EMR
- Highly effective on scarred or fibrotic tissue
- Unique Distal Traction Pin increases snare durability along with top anchoring
- Clear white ergonomic handle for clear marking
- Diameter wire 0.3 mm

Supplier	Product Code	Snare Width (mm)	Description	Catheter		QOM	FAL Code
				Diameter (mm)	Length (cm)		
Vytil	CJ-ADR-23-230-MFO10	10	Hot & Cold Snare	2.3	230	10	
Vytil	CJ-ADR-23-230-MFO15	15	Hot & Cold Snare	2.3	230		
Vytil	CJ-ADR-23-230-MFO20	20	Standard EMR Snare	2.3	230		
Vytil	CJ-ADR2-23-230-O25B	15 / 25	Two in One	2.3	230		FAL23626