

NovaClip R2

Haemostasis



- Innovative handle enabling a 360° rotation with one hand
- Repositionable
- Optimal flexibility of the distal part
- Short tip after release: 7 mm (up to 13 mm max. open) to 9 mm (16 mm max. open)
- MRI compatible

Supplier	Product Code	Maximal Opening (mm)	Jaws Length (mm)	Jaws Angle	Diameter (mm)	Length (cm)	Handle Colour	QOM	FAL Code
Vytil	CJ-HC09-25-230	9	5	135°	2.5	230	■	10	FAL24588
Vytil	CJ-HC11-25-160	11	6	135°	2.5	160	■		FAL24589
Vytil	CJ-HC11-25-230	11	6	135°	2.5	230	■		FAL24590
Vytil	CJ-HC11-25-260	11	6	135°	2.5	260	■		FAL24591
Vytil	CJ-HC13-25-230	13	7	135°	2.5	230	■		FAL23515
Vytil	CJ-HC1690-25-230	16	9	90°	2.5	230	■		FAL23559