


NovaClip R4



- Innovative handle enabling a 360° rotation with one hand
- Repositionable
- Short tip after release: 7 mm (12 clip) and 8 mm (15 clip)
- Optimal flexibility of the distal part
- MRI compatible

Supplier	Product Code	Maximal Opening (mm)	Jaws Angle	Diameter (mm)	Length (cm)	Handle Colour	QOM	FAL Code
Vytil	VTL/HR-230-12-135-Y	12	135°	2.6	230	■	10	FAL23635

Vytil	VTL/HR-230-15-135-Y	15	135°	2.5	230	■	10	FAL23624
-------	---------------------	----	------	-----	-----	---	----	----------