



Stone Extraction Basket - Rotatable and wire-guided

- Nitinol wires with shape memory (4 single-strand)
- Rotary spiral shape
- Compatible with short or long guide wire 0.035" maximum

Supplier	Product Code	Basket diameter (mm)	Catheter Length (cm)	Diameter (mm)	Minimal operating channel diameter (mm)	QOM	FAL Code
Vytil	NLS/BKS-20-24-195G	20	195	3.1	3.7	5	FAL23620
Vytil	NLS/BKS-30-24-195G	30	195	3.1	3.7	5	FAL23725