


**BIÇAKCILAR***Touching Lives One Device at a Time*

# Medical Disposable Solutions

**12<sup>th</sup> Edition**

All products and product specifications identified in the product catalogue are based upon the information available to Bıçakcılar at the time of publication. Bıçakcılar reserves the right to discontinue any of these products or to change any such specifications without prior notice.

## Patients are at the center of every step we take











## Suction Catheter

### Beveled Tip

- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Atraumatic, soft rounded, beveled open tip
- Two lateral eyes with smooth edges
- Color coded connector
- Length: 50 cm (40 and 60 cm are available)



Beveled Tip

| Ref        | Size  | Color   |
|------------|-------|---|
| 190 0601 1 | 06 Ch | Light Green  |
| 190 0801 1 | 08 Ch | Blue         |
| 190 1001 1 | 10 Ch | Black        |
| 190 1201 1 | 12 Ch | White        |
| 190 1401 1 | 14 Ch | Green        |
| 190 1601 1 | 16 Ch | Orange       |
| 190 1801 1 | 18 Ch | Red          |
| 190 2001 1 | 20 Ch | Yellow       |









## Suction Catheter

### Ideal Tip

- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Ideal tip reducing the risk of tissue damage
- Atraumatic, soft rounded, open tip
- Two lateral eyes with smooth edges
- Color coded connector
- Length: 50 cm (40 and 60 cm are available)



Ideal Tip

| Ref        | Size  | Color   |
|------------|-------|---|
| 190 0685 1 | 06 Ch | Light Green  |
| 190 0885 1 | 08 Ch | Blue         |
| 190 1085 1 | 10 Ch | Black        |
| 190 1285 1 | 12 Ch | White        |
| 190 1485 1 | 14 Ch | Green        |
| 190 1685 1 | 16 Ch | Orange       |
| 190 1885 1 | 18 Ch | Red          |
| 190 2085 1 | 20 Ch | Yellow       |

## Suction Catheter

### Ring Tip

Frosted surface for easy insertion

Appropriate degree of hardness provides kink resistance

Open ring tip for bronchial tissue protection

Four small lateral eyes to prevent baro-trauma

Color coded soft connector

Length: 50 cm (40 and 60 cm are available)



Ring Tip

| Ref        | Size  | Color  |
|------------|-------|--------|
| 190 1095 1 | 10 Ch | Black  |
| 190 1295 1 | 12 Ch | White  |
| 190 1495 1 | 14 Ch | Green  |
| 190 1695 1 | 16 Ch | Orange |

## Suction Catheter

### Coudé Tip

Angled tip for endobronchial suction procedures, providing easier access to the right or left bronchus

Frosted surface for easy insertion

Appropriate degree of hardness provides kink resistance

Ideal tip reducing the risk of tissue damage

Atraumatic, soft rounded, open tip

Two small relieving eyes

Marker on funnel connector indicating direction of angled tip

Color coded connector

Length: 50 cm



Coudé Tip

| Ref        | Size  | Color  |
|------------|-------|--------|
| 190 1065 1 | 10 Ch | Black  |
| 190 1265 1 | 12 Ch | White  |
| 190 1465 1 | 14 Ch | Green  |
| 190 1665 1 | 16 Ch | Orange |



Suction Catheter

With Vakon Connector, Ideal Tip

- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Metric type has accurate positioning depth marks at each cm starting 4 cm from the distal end
- Ideal tip reducing the risk of tissue damage
- Atraumatic, soft rounded, open tip
- Two small relieving eyes
- Vakon connector
  - Transparent for easy observation of secretions
  - Color coded
  - Prevents finger contamination
- Length: 50 cm (40 and 60 cm are available)



Ideal Tip

| Metric     | Ref        |       | Size        | Color   |  |
|------------|------------|-------|-------------|---|--|
|            |            |       |             |   |  |
| 019 1190 1 | 190 0551 1 | 05 Ch | Grey        |  |  |
| 019 1191 1 | 190 0651 1 | 06 Ch | Light Green |  |  |
| 019 1192 1 | 190 0751 1 | 07 Ch | Violet      |  |  |
| 019 1193 1 | 190 0851 1 | 08 Ch | Blue        |  |  |
| 019 1194 1 | 190 1051 1 | 10 Ch | Black       |  |  |
| 019 1195 1 | 190 1251 1 | 12 Ch | White       |  |  |
| 019 1196 1 | 190 1451 1 | 14 Ch | Green       |  |  |
| 019 1197 1 | 190 1651 1 | 16 Ch | Orange      |  |  |
| 019 1198 1 | 190 1851 1 | 18 Ch | Red         |  |  |










Suction Catheter

With Vakon II Connector, Ideal Tip

- Frosted surface for easy insertion
- Appropriate degree of hardness provides kink resistance
- Metric type has accurate positioning depth marks at each cm starting 4 cm from the distal end
- Ideal tip reducing the risk of tissue damage
- Atraumatic, soft rounded, open tip
- Two small relieving eyes
- Vakon connector
  - Color coded
  - Prevents finger contamination
- Length: 50 cm (40 and 60 cm are available)



Ideal Tip

| Metric     | Ref        |       | Size        | Color   |  |
|------------|------------|-------|-------------|---|--|
|            |            |       |             |   |  |
| 019 1290 1 | 190 0581 1 | 05 Ch | Grey        |  |  |
| 019 1291 1 | 190 0681 1 | 06 Ch | Light Green |  |  |
| 019 1292 1 | 190 0781 1 | 07 Ch | Violet      |  |  |
| 019 1293 1 | 190 0881 1 | 08 Ch | Blue        |  |  |
| 019 1294 1 | 190 1081 1 | 10 Ch | Black       |  |  |
| 019 1295 1 | 190 1281 1 | 12 Ch | White       |  |  |
| 019 1296 1 | 190 1481 1 | 14 Ch | Green       |  |  |
| 019 1297 1 | 190 1681 1 | 16 Ch | Orange      |  |  |
| 019 1298 1 | 190 1881 1 | 18 Ch | Red         |  |  |



## Suction Catheter

### With Vacuum Control Connector

Frosted surface for easy insertion

Appropriate degree of hardness provides kink resistance

Atraumatic, soft rounded, beveled open tip

Two lateral eyes with smooth edges

Color coded vacuum control connector

Length: 50 cm (40 and 60 cm are available)



Beveled Tip

| Ref        | Size  | Color   |
|------------|-------|---|
| 191 0611 1 | 06 Ch | Light Green  |
| 191 0811 1 | 08 Ch | Blue         |
| 191 1011 1 | 10 Ch | Black        |
| 191 1211 1 | 12 Ch | White        |
| 191 1411 1 | 14 Ch | Green        |
| 191 1611 1 | 16 Ch | Orange       |
| 191 1811 1 | 18 Ch | Red          |

## Suction Catheter

### With Kapkon Connector

Frosted surface for easy insertion

Appropriate degree of hardness provides kink resistance

Atraumatic, soft rounded, beveled open tip

Two lateral eyes with smooth edges

Kapkon connector


Connector cap provides possibility of constant suction

Color coded

Length: 50 cm (40 and 60 cm are available)



Beveled Tip

| Ref        | Size  | Color   |
|------------|-------|---|
| 189 0611 1 | 06 Ch | Light Green  |
| 189 0811 1 | 08 Ch | Blue         |
| 189 1011 1 | 10 Ch | Black        |
| 189 1211 1 | 12 Ch | White        |
| 189 1411 1 | 14 Ch | Green        |
| 189 1611 1 | 16 Ch | Orange       |
| 189 1811 1 | 18 Ch | Red          |

## Kapkon Connector

- Transparent for easy observation of secretions
- Connector cap provides possibility of constant suction
- Accepts funnel connector 4 Ch - 18 Ch



| Ref        |         |
|------------|---------|
| 240 0001 1 | Sterile |
| 240 0001 3 | Bulk    |

## Meconium Aspiration Connector

- Designed to suctioning meconium
- Suction control thumb port for intermittent suctioning
- Attaches to neonatal endotracheal tube connector



| Ref        |  |
|------------|--|
| 236 1001 1 |  |

## Spigot

Closure for catheter funnel connectors  
from 6 Ch to 35 Ch



Spigot I



Spigot II

| Ref        |            |         |
|------------|------------|---------|
| Spigot I   | Spigot II  |         |
| 238 0001 1 | 238 0011 1 | Sterile |
| 238 0001 3 | 238 0011 3 | Bulk    |

## Conical Connector

Luer lock proximal end  
Length: 40 mm



| Ref        |         |
|------------|---------|
| 235 0001 1 | Sterile |
| 235 0001 3 | Bulk    |

Trachea Suction Set

Male/Funnel Connector

- Container for sampling of secretion from tracheal or bronchial tubes
- Transparent container provides clear observation of secretions
- Extra cap ensuring safe sealing of sample container
- Variety of connectors available to fit different suction catheters
- Self-adhesive label for patient identification info



| Ref        |
|------------|
| 187 2001 1 |

Trachea Suction Set

Funnel/Male Connector

- Container for sampling of secretion from tracheal or bronchial tubes
- Transparent container provides clear observation of secretions
- Extra cap ensuring safe sealing of sample container
- Variety of connectors available to fit different suction catheters
- Self-adhesive label for patient identification info



| Ref        |
|------------|
| 187 2101 1 |



## Trachea Suction Set

### Funnel/Kapkon Connector

Container for sampling of secretion from tracheal or bronchial tubes

Transparent container provides clear observation of secretions

Extra cap ensuring safe sealing of sample container

Variety of connectors available to fit different suction catheters

Self-adhesive label for patient identification info



Ref

187 2201 1

## Trachea Suction Set

### Funnel/Kapkon, With Extra Container

Container for sampling of secretion from tracheal or bronchial tubes

Transparent container provides clear observation of secretions

Extra cap ensuring safe sealing of sample container

Variety of connectors available to fit different suction catheters

Self-adhesive label for patient identification info



Ref

187 2301 1