

## Samarit GLIDEBOARD for Seated patient transfer

SAMARIT sliding boards are Nursing aids for moving a seated patient with limited forces from the middle or edge of the bed to a wheelchair, a car seat or a toilet.

Here the patient may not be able to stand but can change position when sitting independently or with the help of a healthcare professional.

While the smooth surface reduces friction resistance to a minimum, non-slip adhesive strips ensure that the GLIDEBOARD doesn't slip away. The GLIDEBOARD is made of an unbreakable material, is acid-resistant and easy to clean and disinfect. Five versions are available for different purposes.

The SAMARIT GLIDEBOARD range has been developed for special use with the support of therapists.

### Glideboard - THE GREAT



Model-No.:	476.0080
UDI/GTIN:	7640144460938
Use:	to transfer patients from the edge of the bed to a wheelchair. The product is used mainly today to transfer patients from wheelchairs into large vehicles (SUVs etc.)
Use by:	Gap < 30 cm / 11 inch = 130 kg Gap < 20 cm / 9 inch = 150 kg
Size:	80 x 33 cm
Weight:	2.6 kg
Use by:	Gap < 30 cm / 11 inch = 130 kg Gap < 20 cm / 9 inch = 150 kg
Weight Limit:	220 kg

### Glideboard - CURVED “Banana”



Model-No.:	476.0012
UDI/GTIN:	7640144460945
Use:	to enable sedentary patients to transfer independently to wheelchairs, car seats or chairs.
Size:	68 x 23 cm
Weight:	1.2 kg
Weight Limit:	200 kg