


The perfect balance of strenght & flexibility



- 21 Gauge Introducer Needle with Echogenic Tip
 - Super Sharp Needle
- Guidewires
 - Nitinol Material
- Tearaway Sheet
 - Wide mouth opening

Reference	Description	Order Qty
MRCTP41028	4F x 55cm PRO-PICC®, Single Lumen	5pcs/box
MRCTP51028	5F x 60cm PRO-PICC®, Single Lumen	5pcs/box
MRCTP52028	5F x 55cm PRO-PICC®, Double Lumen	5pcs/box

Set contains:
 1 PICC w/ Stylet
 1 Peelable Sheath Introducer
 1 Introducer Noodle
 1 70cm Guidewire
 1 Mini Scalpel
 1 Tape Measure
 1 Catheter Securement Device

Medcomp is distributed by:

Duomed Belgium • Oeyvaersbosch 12 • 2630 Aartselaar • Belgium
 +32 3 870 11 11 • info@duomed.com • www.duomed.com