

ProFUSE^{CT}

Power injectable CT implantable port



RAISED RIM

The raised rim provides a clear tactile response during palpation to assist in septum location for infusion therapy.



LOCKING MECHANISMS

Clear locking mechanisms provide multiple sensory confirmations, through audible, tactile, and visual responses. All three features afford added security for the healthcare professional to ensure proper securement of the catheter to the port.

CT IDENTIFICATION

Utilizing a unique printing technique, "CT" can be visualized under x-ray, to provide recognition of the power injectability of the Pro-Fuse CT device. All patients with the Pro-Fuse^{CT} receive a patient packet, which includes additional identification tools.



Reference	Description	Order Qty
MRCTT80001	Profuse Standard	5pcs/box

Set contains:

- | | | |
|--------------------------|--|----------------------------------|
| 1 Catheter | 1 .035" x 70cm "J" Marked Guidewire w/ Advancer | 1 Blunt Tip Needle |
| 1 CT Implantable Port | 1 10cc Syringe | 1 Valved Tearaway Introducer IFU |
| 2 Catheter Locks | 1 Valved Tearaway Introducer | 1 Port IFU |
| 1 Scalpel | 1 Tunneler | 1 Patient Information Pack |
| 1 18ga Introducer Needle | 2 22ga Huber Needles (1-Straight, 1-Right Angle) | 1 Patient Chart Sticker |
| 1 Vein Pick | | |

Medcomp is distributed by:

Duomed Belgium • Oeyvaersbosch 12 • 2630 Aartselaar • Belgium
 +32 3 870 11 11 • info@duomed.com • www.duomed.com