




SonoScape

DEFINED MOVEMENT, UNLIMITED PURPOSE

EC-P560 Series,
SonoScape's Slim Colonoscopes with ATS and PFT Technology





The versatility of SonoScape's slim scopes originates from their well-balanced configurations.

With a small distal end, a slim body, and a satisfying instrument channel size, not only do they excel in routine diagnostic and therapeutic procedures, but are the go-to scopes for facing colon narrowing, stenosis and volvulus, or geriatric patients etc. In addition, their excellent field and depth of view assures the high quality of a colonoscopy.

Equipped with waterjet, ATS* and PFT† technology, SonoScape's slim scopes are physician focused. Their stiffness can be adjusted in continuous levels, and to the perfect nuance. The scope tubing is improved so it becomes more sensitive to the endoscopist's hand gestures, even easier to manoeuvre under challenging situations.

*ATS: Alterable Tubing Stiffness

†PFT: Precise Force Transmission

EC-P560 Series

	EC-P560	EC-P560T
Field of view (°)	170	
Depth of view (mm)	2-100	
Distal-end diameter (mm)	10.8	
Insertion tube diameter (mm)	11.5	
Instrument channel diameter (mm)	3.2	
Angulation (°)	Up/Down 180, Left/Right 160	
Working length (mm)	1350	1700
Waterjet	Available	
Special features	ATS (Alterable Tubing Stiffness), PFT (Precise Force Transmission)	
Connection	OneTouch, waterproof	
Compatibility	HD-550	

Exclusive distributor of SonoScape in The Netherlands



Sales
+31 318 54 32 23
sales.nl@duomed.com

Service
+31 318 51 89 48
service.nl@duomed.com

The team

