




**SonoScape**

# VIDEO GASTROSCOPE

## EG-550/EG-550L



# Boosted Maneuverability



## Ergonomic experience

A new scope tubing paired with a light-weight handle, the upgraded 550 series scopes offer great comfort and ease during routine and complex operations, helping to shorten the procedure time, and to reduce both patients' discomfort and yours.

## Fulfilling visual field

After an easy intubation, observe to full extent owing to the scope's short bending radius and large field of vision. Diagnose and treat with confidence, thanks to its water jet function and generous instrument channel.



## Ingenious design

One-Touch connection, and real water-proof, the 550 series scopes can be plugged and unplugged, cleaned and disinfected straightforwardly at a minimal risk of water damage by mismanipulation, which indeed lightens your team's workload.

### 550 SERIES SCOPES

	EG-550	EG-550L
Field of view	140°	140°
Depth of view	3-100mm	3-100mm
Distal End Dia.	9.3mm	9.8mm
Insertion Dia.	9.3mm	9.8mm
Instrument Dia.	2.8mm	3.2mm
Angulation	U210°, D90°, L/R100°	
Working length	1050mm	1050mm
Water Jet	Yes	Yes
OneTouch	Yes	Yes
Optical Fiber Trans	Yes	Yes

Exclusive distributor of SonoScape in The Netherlands



**Sales**  
+31 318 54 32 23  
sales.nl@duomed.com

**Service**  
+31 318 51 89 48  
service.nl@duomed.com

The team

