



SonoScape



EG-550Z

Optical Magnification Gastroscope





A distal end diameter as slim as 5.4 mm

makes it easier for both Healthcare Professionals and patients



Waterproof enclosure design



140°

Wide field of view of 140° dual light guides

resulting in a brighter and better endoscopic observation



One-Step connection



Large instrument channel of 2.2 mm

to achieve an adequate suction capacity

EG-P550N

Field of view	140°
Depth of field	3 – 100 mm
Angulation	Up 210° / Down 90°
	Right 100° / Left 100°
Distal-end outer diameter	5.4 mm
Insertion tube outer diameter	5.8 mm
Instrument channel inner diameter	2.2 mm
Working length	1100 mm

Exclusive distributor of SonoScape in The Netherlands



Sales
+31 318 54 32 23
sales.nl@duomed.com

Service
+31 318 51 89 48
service.nl@duomed.com

The team

