

Complete Drainage System

SKATER™



DUOMED




ARGON
MEDICAL DEVICES

Formulated to promote r

Large skived holes
help maximise drainage

Position markers
confirm catheter position
after placement


**Locking
shoulder tip**
reduces buckling
during insertion

Removable pigtail s

SLIP-COAT[™] hydrophilic coating
significantly reduces friction for
easier catheter insertion



maximum drainage flow



straightener

Choice Lock™
trocar stylet
facilitates one-handed
placement

French size indicator
is clearly marked

**Consistent catheter-to-hub
inner lumen size**
prevents
bottlenecking

**Low profile
locking system
with snap-off tab**
provides patient
comfort

The SKATER™ Drainage System provides a complete portfolio of catheters and accessories for all purpose, nephrostomy, biliary and centesis procedures.

All Purpose/Nephrostomy Drainage

Catheters for percutaneous drainage of abscesses and fluid collections using the direct stick technique or kidney using the Seldinger technique

- Expansive size offering of 6F-16F and 15cm-60cm
- Available in locking and non-locking configurations
- Choice Lock™ trocar stylet allows for single-handed abscess drainage placement (not available on all sizes)
- Designed for Seldinger technique
- Can be used safely with ethyl alcohol and doxycycline

Also Available:

- SKATER Introducer Sets
- SKATER Kits
- SKATER Centesis
- SKATER Mini-Loop
- Introducer Sheaths/Needles (ISN)
- Hawkins™ Blunt Needle Access System
- Drainage Bags
- SKATER Fix Catheter Fixation Device

Biliary Drainage

Catheters for drainage of the biliary tract

- Radiopaque marker band provides clear view of the most proximal drainage hole for accurate placement in the biliary tract
- Additional drainage holes for biliary duct to duodenum drainage
- Can be used safely with ethyl alcohol and doxycycline



Ordering Information

Each Set Includes:

- 1 x SKATER Catheter
- 1 x Removable pigtail straightener
- 1 x Metal Stiffening Cannula
- 1 x Flexible Stiffening Cannula (Not available on 6F and 7F)
- 1 x Locking system with snap-off tab (not available on non-locking pigtail catheters)

SKATER All Purpose and Nephrostomy Drainage Set Locking (Box of 5)				
CATALOGUE NO.	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET	TROCAR SIZE
756606015	6F x 15cm	.035"	Yes	19G
756606020	6F x 20cm	.035"	Yes	19G
756606025	6F x 25cm	.035"	Yes	19G
756606035	6F x 35cm	.035"	No	
756607015	7F x 15cm	.035"	Yes	18G
756607020	7F x 20cm	.035"	Yes	18G
756607025	7F x 25cm	.035"	Yes	18G
756608015	8F x 15cm	.038"	Yes	17G
756608020	8F x 20cm	.038"	Yes	17G
756608025	8F x 25cm	.038"	Yes	17G
756608030	8F x 30cm	.038"	Yes	17G
756608035	8F x 35cm	.038"	No	
756608045	8F x 45cm	.038"	No	
756610015	10F x 15cm	.038"	Yes	17G
756610020	10F x 20cm	.038"	Yes	17G
756610025	10F x 25cm	.038"	Yes	17G
756610030	10F x 30cm	.038"	Yes	17G
756610035	10F x 35cm	.038"	No	
756610045	10F x 45cm	.038"	No	
756610060	10F x 60cm	.038"	No	
756612015	12F x 15cm	.038"	Yes	14G
756612020	12F x 20cm	.038"	Yes	14G
756612025	12F x 25cm	.038"	Yes	14G
756612030	12F x 30cm	.038"	Yes	14G
756612035	12F x 35cm	.038"	No	
756612045	12F x 45cm	.038"	No	
756612060	12F x 60cm	.038"	No	
756614015	14F x 15cm	.038"	Yes	14G
756614020	14F x 20cm	.038"	Yes	14G
756614025	14F x 25cm	.038"	Yes	14G
756614030	14F x 30cm	.038"	Yes	14G
756614035	14F x 35cm	.038"	No	
756614045	14F x 45cm	.038"	No	
756614060	14F x 60cm	.038"	No	
756616025	16F x 25cm	.038"	Yes	14G
756616030	16F x 30cm	.038"	Yes	14G
756616035	16F x 35cm	.038"	No	
756616045	16F x 45cm	.038"	No	
756616060	16F x 60cm	.038"	No	

*Flexible stiffening cannula not available on 6F and 7F.

*Additional ordering information is available on the back.

SKATER All Purpose and Nephrostomy Drainage Set Non-Locking (Box of 5)				
CATALOGUE NO.	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET	TROCAR SIZE
756106015	6F x 15cm	.035"	Yes	19G
756106020	6F x 20cm	.035"	Yes	19G
756106025	6F x 25cm	.035"	Yes	19G
756107020	7F x 20cm	.035"	Yes	18G
756107025	7F x 25cm	.035"	Yes	18G
756108020	8F x 20cm	.038"	Yes	17G
756108025	8F x 25cm	.038"	Yes	17G
756108030	8F x 30cm	.038"	Yes	17G
756110020	10F x 20cm	.038"	Yes	17G
756110025	10F x 25cm	.038"	Yes	17G
756110030	10F x 30cm	.038"	Yes	17G
756112020	12F x 20cm	.038"	Yes	14G
756112025	12F x 25cm	.038"	Yes	14G
756114020	14F x 20cm	.038"	Yes	14G
756114025	14F x 25cm	.038"	Yes	14G
756116025	16F x 25cm	.038"	Yes	14G

*Flexible stiffening cannula not available on 6F and 7F.

SKATER Mini-Loop Drainage Catheter w/Locking Pigtail (Box of 5)				
CATALOGUE NO.	FRENCH x LENGTH	MAX GUIDEWIRE SIZE	TROCAR STYLET	TROCAR SIZE
758606015	6F x 15cm	.035"	Yes	19G
758606020	6F x 20cm	.035"	Yes	19G
758606025	6F x 25cm	.035"	Yes	19G
758607015	7F x 15cm	.035"	Yes	18G
758607025	7F x 25cm	.035"	Yes	18G
758608015	8F x 15cm	.035"	Yes	18G
758608025	8F x 25cm	.035"	Yes	18G
758610015	10F x 15cm	.038"	Yes	17G
758610025	10F x 25cm	.038"	Yes	17G
758612015	12F x 15cm	.038"	Yes	14G
758612025	12F x 25cm	.038"	Yes	14G
758614025	14F x 25cm	.038"	Yes	14G

SKATER Biliary Drainage Catheter w/Locking Pigtail (Box of 5)		
CATALOGUE NO.	FRENCH x LENGTH	MAX GUIDEWIRE SIZE
755308040	8F x 40cm	.038"
755310040	10F x 40cm	.038"
755312040	12F x 40cm	.038"

SKATER Biliary Drainage Catheter w/Non-Locking Pigtail (Box of 5)		
CATALOGUE NO.	FRENCH x LENGTH	MAX GUIDEWIRE SIZE
755208040	8F x 40cm	.038"
755210040	10F x 40cm	.038"
755212040	12F x 40cm	.038"



DUOMED

Duomed Swiss AG
 Grenzstrasse 5a, 6214 Schenkon
 Tel. +41 (0)41 510 07 00
www.duomed.com



© 2021 Argon Medical Devices, Inc. All Rights Reserved.
 No reproduction of any portion of this work shall be made without prior permission of the author.



96-6300-24/0521C