

Duomed Swiss AG

Grenzstrasse 5a, 6214 Schenkon

Tel. +41 (0)41 510 07 00

www.duomed.com




 Claritron®

7.5Fr

Single-use Ureteroscope

7.5Fr SUV-2C-B Single-use Ureteroscope



Superfine Shaft More Efficient Irrigation



Irrigation Up By 150 % (Compared with 9.0Fr Ureteroscope in 12Fr Ureteral Access Sheath)

Thinner and uniform 7.5Fr access shaft, which allows for increased flow of irrigation fluid around the ureteroscope to maintain lower intrarenal pressure.



Easier Insertion

The 7.2Fr distal end and 7.5Fr shaft are designed for easier insertion and less traumatic for patients. Also, it is more suitable for use with children.



High Performance

| Additional $\pm 60^\circ$ Rotation Function
| Excellent Passive Bending

Product Code	SUV-2C-B
Outer Diameter of Insertion Section	7.5Fr/2.5mm
Outer Diameter of distal End	7.2Fr/2.4mm
Working Length	680mm
Total Length	990mm
Working Channel (ID)	3.6 Fr/1.2mm
Field of View	120°
Depth of Field	0.5-120mm
Deflection	285° Up/Down



SCIVITA MEDICAL TECHNOLOGY CO.,LTD.

+86-512-8187 7788

www.scivitamedical.com

info@scivitamedical.com

No. 2, Qingqiu Street, Suzhou Industrial Park, 21500
0 Suzhou, Jiangsu Prov., P.R.CHINA

