

SmartBand[®]

Multi-Band Ligation System



SmartBand[®]

Multi-Band Ligation System



Specifically designed to improve ligation performance...

- ▶ The SmartBand® Multi-Band Ligation System is used to endoscopically ligate esophageal varices at or above the gastroesophageal junction and to ligate internal hemorrhoids.


The SmartBand Multi-Band Ligation Kit...

Consists of a ligation handle with universal connector, loading device, flush tube, Pentax adaptor, and barrel with 6 bands.


The bands are...

- Designed to deliver maximum tissue compression and gripping force
- Uniquely manufactured and packaged to enable a 24 month shelf life


SmartBand Multi-Band Ligation Kit



SmartBand SafeGrip Multi-Band Ligation Kit



SmartBand Ligation bands delivered higher average compression forces than competition during third party testing¹



SmartBand SafeGrip Ligation bands delivered 33% higher gripping force than competition during third party testing¹

The SmartBand Multi-Band Ligation System...

Features a universal connector that attaches the handle to the endoscope and is compatible with multiple endoscope styles




Includes 6 backloaded bands on a clear, smooth barrel for optimal visualization



Offers a 10FR inner diameter handle stem that enables maximum suction capabilities



Provides a dedicated flush port for clearing the working channel during the procedure



The SmartBand Packs offer a reloadable barrel of bands for procedures that require more than 6 bands – helping to reduce costs



SmartBand[®]
Multi-Band Ligation System

For more information on the SmartBand Multi-Band Ligation System, contact STERIS Endoscopy today or visit www.steris.com

Product Number	Description	Endoscope Diameter Compatibility (mm)	Unit of Measure
SLK6*	SmartBand Multi-Band Ligation Kit - Components: Ligation Handle with Universal Connector, Loading device, Flush Tube, Pentax Adaptor, and Barrel with 6 Bands	8.6 - 11.6	EA
SLK6LF**	SmartBand SafeGrip Multi-Band Ligation Kit - Components: Ligation Handle with Universal Connector, Loading device, Flush Tube, Pentax Adaptor, and Barrel with 6 Bands	8.6 - 11.6	EA
SLP6*	SmartBand Multi-Band Ligation Pack - Components: Deployment Cord with Barrel of 6 Bands	8.6 - 11.6	EA
SLP6LF**	SmartBand SafeGrip Multi-Band Ligation Pack - Components: Deployment Cord with Barrel of 6 Bands	8.6 - 11.6	EA

*This product is made with natural rubber latex.

**This product is not made with natural rubber latex.

This device is supplied non-sterile and is disposable, single-use only.



DUOMED

Duomed Swiss AG

Grenzstrasse 5a, 6214 Schenkon

Tel. +41 (0)41 510 07 00

www.duomed.com



SCAN to view
full set-up video.



[1] Third party testing data on file

© 2021 STERIS Endoscopy or its affiliates. All rights reserved.
Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by STERIS or its affiliates.

Follow us:    



STERIS

Endoscopy

STERIS

5960 Heisley Road

Mentor, OH 44060-1834 – USA

Customer Service: 800-548-4873

Fax: 440-639-4495

www.steris.com