

JAWZTM

ENDOMYOCARDIAL BIOPSY FORCEPS




CUTTING EDGE CARE

Precise | Sharp | Safe



CUTTING EDGE CARE

Precise | Sharp | Safe



**Precision-molded
single-use design**

To provide a solid
"instrument" feel

Ergonomic handle

Maintains desired thumb
position while providing
comfort

Enhanced control

With full rotation

**Advanced
dual cutting head**

Designed for sharper,
cleaner cuts

**Balance of rigidity
and flexibility**

To advance into the tissue
and safely navigate
to the target


Improved Tissue Access

Tailored to your technique,
JAWZ simplifies tissue
retrieval with formable and
non-formable options in
stainless steel or FEP-coated
finishes. Choose from
Straight, Pre-Curved or
Maxi-Curved configurations

¹ Jawz Endomyocardial Biopsy Forceps from Argon Medical: A Review of Safety and Efficacy Based on Literature Review and Bench Model Comparison With the Cordis Biopptome. P. Stibbs, J. Reinault, M. Parsons

Precise, Sharp, Safe

JAWZ™ offers a clear-cut approach to better endomyocardial biopsy samples. It is designed to provide higher quality samples with fewer samples required. The portfolio consists of disposables that maintain sharp dual cutting heads with design flexibility.



PRECISE

SHARP

SAFE

AT HEART

Allowing superior
surface alignment

Ensuring diagnostic quality,
higher tissue volume and
less crush artifact

Low complication rates¹

Straight

Box of 1



Catalogue #	Jaw Size	Jaw OD (French)	Jaw Volume	Length	Curve	Coating	Formability
190090	1.5mm	5F	1.80mm ³	50cm	Straight	Stainless Steel	Formable
190095	1.5mm	5F	1.80mm ³	105cm	Straight	Stainless Steel	Formable
190060	1.8mm	6F	2.46mm ³	50cm	Straight	Stainless Steel	Formable
190065	1.8mm	6F	2.46mm ³	50cm	Straight	Stainless Steel	Non-Formable
190066	1.8mm	6F	2.46mm ³	70cm	Straight	FEP	Formable
190070	1.8mm	6F	2.46mm ³	105cm	Straight	Stainless Steel	Formable
190075	1.8mm	6F	2.46mm ³	105cm	Straight	Stainless Steel	Non-Formable
190085	2.2mm	7F	4.97mm ³	105cm	Straight	Stainless Steel	Formable
190086	2.2mm	7F	4.97mm ³	105cm	Straight	FEP	Formable
190020	2.4mm	7.5F	5.43mm ³	50cm	Straight	Stainless Steel	Formable
190040	2.4mm	7.5F	5.43mm ³	105cm	Straight	Stainless Steel	Formable

Pre-curved

Box of 1



Catalogue #	Jaw Size	Jaw OD (French)	Jaw Volume	Length	Curve	Coating	Formability
190080	2.2mm	7F	4.97mm ³	50cm	Pre-Curved	FEP	Formable
190030	2.4mm	7.5F	5.43mm ³	50cm	Pre-Curved	FEP	Formable

Maxi-Curved

Box of 1



Catalogue #	Jaw Size	Jaw OD (French)	Jaw Volume	Length	Curve	Coating	Formability
190051	1.8mm	6F	2.46mm ³	50cm	Maxi-Curved	FEP	Formable
190081	2.2mm	7F	4.97mm ³	50cm	Maxi-Curved	FEP	Formable
190031	2.4mm	7.5F	5.43mm ³	50cm	Maxi-Curved	FEP	Formable

Explore Cutting Edge Care



Indication for Use: The 50cm and 70cm disposable radiopaque endomyocardial biopsy forceps are designed for right ventricular biopsies using the jugular approach, and for right or left ventricular pediatric biopsies using the femoral approach. The 105cm disposable radiopaque endomyocardial biopsy forceps are designed for right or left ventricular biopsies using the femoral approach.

Contraindications: There are no medical conditions (contraindications) associated with the Endomyocardial Biopsy Forceps. Assessment of patient overall physical state by the physician is required prior to use of the device. **RX only**

ArgonMedical.com

© 2024 Argon Medical Devices, Inc. All rights reserved. Argon Medical Devices, the stylized logo, and the product and service marks mentioned herein are trademarks or registered trademarks of Argon Medical Devices, Inc. in the United States and other countries.

ARGON
MEDICAL DEVICES