

Safety in a new light



SecuFix Uterus Manipulator


Pioneering System with
Colpotransillumination and Vacuum Fixation

SecuFix Uterus Manipulator

Safety in a new light

The development of patient-friendly, minimally invasive hysterectomy techniques and the associated good results, now mean that hysterectomy can today be carried out laparoscopically in many cases.

The SecuFix Uterus Manipulator defines new benchmarks for safety. The atraumatic fixation of the uterus using vacuum aspiration and "Colpotransillumination", i.e. the option of visualizing tissue using distal, annular illumination for the cervical cap are the significant new product features. They are important factors for an optimum surgical result and permit compliance with oncological criteria, e.g. the avoidance of intraperitoneal tumor cell dislocation.



SecuFix Uterus Manipulator Starter Kit

SecuFix Uterus Manipulator	8370701
SecuFix Uterus Probe with longitudinal scale	8370702
SecuFix Colpotransilluminator Vacuum Tube, steril (pack of 10 sets, packed singly, sterile)	41708370
Sealing Cap (pack = 10 pcs).....	89.02
ENDOLIGHT Flex LED Light Source bndl High Power LED, LED service life > 30,000 h, color temperature approx. 6,500 K, continuously variable electronic light quantity control, VESA 100 interface, 100–240 V a.c., 50/60Hz, Dimensions (W x H x D): 150 x 55 x 220 mm Comprising: ENDOLIGHT Flex LED (5166001) and power cable (244003)	51660011
Fiber Light Cable bndl Fiber bundle Ø 3,5 mm (color marking orange), length 2,3 m	806635231

Accessories

Cleaning Brush	7970705
Cleaning Brush	7970409
RIWO-System-Tray Instrument container, inside dimensions (L x W x H): 500 x 200 x 60 mm, steam sterilizable	33024
Irrigation Tube for machine reprocessing, length 600 mm.....	163915

