



INJECTION NEEDLE

ELITESHARP Injection needle provides cystoscopic injection into urethra, bladder neck and bladder wall.



Features:

- 304 Stainless Steel needle
- Bevel tip for easy penetration
- Ultrasharp Needle Tip
- Injection Depth marking for safe injection
- Luer Lock connector to fit syringes
- Side wing for easy manipulation

PRODUCT NAME	CATALOG NO	DIAMETER (GAUGE)	LENGTH (mm)
Injection Needle, Metal	MIN1835M	18	35
Injection Needle, Metal	MIN1935M	19	35
Injection Needle, Metal	MIN2035M	20	35
Injection Needle, Metal, Plastic Hub	MIN1835	18	35
Injection Needle, Metal, Plastic Hub	MIN1935	19	35
Injection Needle, Metal, Plastic Hub	MIN2035	20	35
Injection Needle, Metal, Plastic, Tapered Tip	MINT1935	19	35
Injection Needle, Plastic	MIN1835P	18	35
Injection Needle, Plastic	MIN1935P	19	35
Injection Needle, Plastic	MIN2035P	20	35

Distributed in Belgium & Luxembourg by:

