

COBRA CLIP

ENDOSCOPIC HEMOCLIPS

SECURE / GRIP / RESILIENCE



One-piece clip with elastic structure

Strong closure force adapted to target tissues

Anti-slip jaw tip design

Excellent grasping during tissue mobilization

Ergonomical handle

Precise opening and closing control for efficient release management. Safety lock engages after release

1:1 rotation


Intuitive handling while keeping the hand on the handle

Flexible catheter end

Reliable closing and rotation performance, even under challenging conditions

Integral color code

Quick identification of opening width in emergency situations



References	Clip		Catheter		Color code
	Maximal opening (mm)	Jaws angle (°)	Diameter (mm)	Length (cm)	
VTL/HC09-135-23-235	9	135	2.3	235	
VTL/HC11-135-23-235	11	135	2.3	235	
VTL/HC11-135-23-165*				165	
VTL/HC13-135-23-235	13	135	2.3	235	
VTL/HC16-90-23-235	16	90	2.3	235	
VTL/HC16-90-23-165*				165	

* Special Request - delivery time on demand



Vytil

161 avenue Gallieni, 93170 Bagnolet, France
info@vytil.com
www.vytil.com