

IOB Convective Warming System

505 warming unit product specifications



Dimensions: H: 29.5 cm x 22 cm x 22 cm

Weight: 5.4kg

Relative Noise Level: ≤52 dB(A)

Operating Temperature:

- High : 43° ± 3°C
- Mid : 38° ± 3°C
- Low : 32° ± 3°C

Device Rating:

- 110-120 VAC, 60 Hz, 10 A
- 220-240 VAC, 50 Hz, 8 A

Filter: HEPA 0.2 µm. Recommended use: Every 12 months or 500 hours

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-1-8, IEC 60601-2-35

Reference: IOB-505



Single-use convective warming blankets

Primary packaging: Pouch

Secondary packaging: Cardboard


Minimum Delivery Quantity (MDQ): 40

PVC & Latex free, components PP and LDPE




Standard blankets

Upper Body
Ref: IOB-003
202cm x 64cm




Multi-Position Upper Body
Ref: IOB-030
198cm x 80cm




Lower Body
Ref: IOB-002
142cm x 120cm




Dual Port Torso
Ref: IOB-019
142cm x 120cm




Full Body
Ref: IOB-004
210cm x 120cm



Multi-Access
Ref: IOB-018
210cm x 120cm




Cath Lab
Ref: IOB-015
180cm x 74cm



Pediatric blankets


Pediatric Full Body

Ref: IOB-007
170cm x 100cm



Large Pediatric Under Body

Ref: IOB-009
160cm x 80cm



Pediatric Lower Body

Ref: IOB-012
64cm x 84cm




Infant Full Body & Under Body

Ref: IOB-025
120cm x 80cm



Infant Under Body


Ref: IOB-028
100cm x 100cm



Specialty blankets

Surgical Access

Ref: IOB-016
210cm x 120cm




Sterile Cardiac

Ref: IOB-020S
142cm x 120cm




Adult Under Body

Ref: IOB-006
215cm x 100cm



Lithotomy & Orthopedic Underbody

Ref: IOB-029
200cm x 100cm



Distributed by

Duomed Scandinavia • Lejrvej 25 • 3500 Værløse • Denmark
T +45 44 47 47 48 • info.dk@duomed.com • www.duomed.com

Imported by

